


Short Instructions

Eppendorf ThermoMixer F0.5/F1.5/F2.0/FP English (EN)

1 Operating controls



- **1 Temperature keys with control LEDs** Select a temperature or switch off temperature control.
- 2 Key short Short Mix runs as long as the short key is pressed.
- **3 Key start/stop** Start or stop mixing/temperature control.
- 4 Arrow keys speed

Set the mixing frequency. Keep the arrow key pressed: quick setting.

5 Arrow keys temp

Set the temperature. Keep the arrow key pressed: quick setting. As soon as the set temperature is changed, the device begins to perform temperature control.

- 6 Display
- 7 Menu arrow key Navigate in the menu: set key lock or volume.
- 8 Key menu/enter Open the menu. Confirm your selection.

Short Instructions

Eppendorf ThermoMixer F0.5/F1.5/F2.0/FP English (EN)

1 2 3 4 (5 BB) 创 U 30/23300 °Ĉ rpm 5 7 6

Display ThermoMixer F0.5/F1.5/F2.0/FP

1 **Device status**

2

C Device is performing mixing/temperature control.

2 ThermoTop

ThermoTop has been attached. In order to 5 prevent the formation of any condensate in a reliable manner, the device first heats the ThermoTop until it reaches the set temperature. The tempering of the thermoblock occurs with a delay. After inserting samples into a pre-heated thermoblock, the displayed actual temperature may fall temporarily.

3 Speaker

• Speaker switched on. **K** Speaker switched off.

4 Key lock

Key lock activated: parameters cannot be changed.

No key lock.

Mixing frequency

6 Actual temperature

When the actual temperature flashes on the display, the device does not operate.

7 Set the temperature

When the set temperature has been reached, only one value is displayed.