## **16510, 16511**

Directions for Polycarbonate  
Filter Holders

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**English – page 2**

Before reading this instruction manual, please fold out this page and the last pages showing the illustrations/photos.

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
## Contents

1. Technical Data
2. How to Unpack
3. Cleaning the Filter Holder
4. Clarifying Filtration
5. Sterilizing Filtration
6. Microbiological Testing
7. Sterility Testing
8. Index to Recommended Accessories

### 1. Technical Data

Material	Makrolon® (polycarbonate) Polypropylene filter support
Sealing	Silicone O-ring, 40 × 5 mm (6980110)
Filter diameter	47 mm
Prefilter diameter	37 mm
Filtration area	12.5 cm <sup>2</sup>
Capacity	250 ml
Connections	Female Luer slip (3 each in the lid of the funnel and receiver flask, 2 on opposite sides of the filter rest [table]). 2 hose nipples 6980234 for 10 mm hose. Outlet of SM 16511: 20 mm outer diameter.
Operating pressure	Vacuum or 2 bar pressure max.
Sterilization	Autoclave at 121 °C

### 2. How to Unpack

Remove filter holder from packaging. Make sure to remove all parts.

### 3. Cleaning the Filter Holder

- 3.1 Completely disassemble filter holder (Fig. 1).
- 3.2 Clean all parts with a laboratory detergent which is commercially available (for metal, glass, and synthetics).
- 3.3 Rinse parts with hot water first and then with distilled water.
- 3.4 Dry parts using compressed air or drip-dry. Use of towels or cloth to dry is not recommended since they may leave lint particles on filter holder. If instructions for cleaning are properly followed, the material of the filter holder will remain clear, and your filter holder will always be ready to use.

#### 4. Clarifying Filtration

- 4.1 Screw lid onto the receiver flask (Fig. 2).
- 4.2 Close off the three connections (female Luer slip) on the lid with the plugs provided (Fig. 3).
- 4.3 Screw filter rest (table) onto the lid of the receiver flask (Fig. 4).
- 4.4 Place the filter support on the filter rest and then a 47 mm  $\varnothing$  membrane filter on top of it (If necessary, place 37 mm  $\varnothing$  prefilter on the membrane filter.) (Fig. 5).
- 4.5 Next, attach the upper portion of the filter holder to the filter rest (Fig. 6).
- 4.6 Screw second lid onto the upper portion (after inserting the proper O-ring) (Fig. 7).
- 4.7 For vacuum filtration, screw the filling funnel onto the upper lid, and insert a plug in one of the lateral connections of the filter rest (Fig. 8).
- 4.8 Connect a vacuum hose to a hose nipple and attach the end of the hose to a vacuum source (water jet pump or a water filter vacuum pump) before inserting the hose nipple in the other lateral connection (Fig. 9)
- 4.9 Pour the liquid to be filter into the filling funnel, and apply vacuum to the bottom portion of the filter holder to filter (Fig. 10).

As an **alternative**, you may also filter using pressure. Follow the same steps as described in 4. -4.6. Both lateral connections on the filter rest are left open.

After pouring the liquid to be filtered into the filling funnel, screw on closing cap and insert a plug into each of the two upper connections. Attach a hose nipple to the third connection (after connecting it first to a pressure hose, and attaching hose to a pressure source). Apply pressure (2 bar max.) to the filter holder in order to filter (Fig. 11).

#### 5. Sterilizing Filtration

- 5.1 Insert a 0.2  $\mu\text{m}$  membrane filter (and prefilter, if required) into the filter holder, and cap all of the connections on the two lids (top and bottom). Screw closing cap onto the filling port. Insert a plug into one of the lateral connections of the filter rest. Attach a syringe filter holder, SM 16514, to the other connection, after previously inserting a completely hydrophobic filter SM 13307 - 013 into the syringe filter holder (Fig. 12).

5.2 Wrap the second closing cap in aluminum foil, and autoclave the assembled filter holder as well as the wrapped closing cap for 30 minutes at 121 °C.

5.3 Perform sterilizing filtration as you would pressure filtration. Observe the usual precautions to avoid contamination.

5.4 After completing filtration, shut off pressure source, and remove the nipple along with the attached pressure hose; simply pull off nipple from the connection (Fig. 13).

5.5 Unscrew the entire upper portion of the filter holder and remove (Fig. 14).

5.6 Cap the receiver flask with the other autoclaved closing cap (Fig. 15).

## **6. Microbiological Testing**

6.1 Insert the appropriate membrane filter, 47 mm Ø, into the filter holder.

6.2 Plug all connections on the top and bottom lids using the plugs provided.

6.3 Insert a plug into one of the lateral connections, and cover the other with aluminum foil.

6.4 Screw filling funnel onto the filling port and cover with aluminum foil.

6.5 Now the filter holder is ready to be autoclaved at 121 °C for 30 min. (Fig. 16).

6.6 After autoclaving, remove the aluminum foil. Connect the filter holder to a vacuum source by attaching a hose nipple to one of the lateral connections and connecting a hose to the nipple and vacuum source.

6.7 Pour the liquid to be filtered into the filling funnel and apply vacuum to the filter holder to begin filtration (Fig. 17).

6.8 After completing filtration, pull the hose nipple along with the hose off the connection, and remove the upper portion of the filter holder (Fig. 18).

6.9 Remove the membrane filter from the filter rest using sterile forceps. Incubate it on a nutrient pad, on agar, or in other culture media (Fig. 19).

## **7. Sterility Testing**

7.1 You may attach the upper portion only of SM 16511 to a vacuum manifold (SM 16826, 3-place), and perform the sterility test (Fig. 20).

## 8. Index to Recommended Accessoires

SM 16623	Rubber vacuum hose (1 m)
SM 16611	Water let pump with $\frac{3}{4}$ " female thread
SM 16612	Laboratory vacuum pump, 220V, 50 Hz
SM 16615	Laboratory vacuum pump, 110 V, 50 Hz
SM 16610	Woulff's bottle, 500 ml
SM 16624	Stainless steel forceps
SM 16826	Stainless steel 3-place manifold
SM 16827	Stainless steel flushing tube for SM 16826, with taps
SM 16829	Stainless steel 6-place manifold
SM 16830	Stainless steel flushing tube for SM 16829, with taps
SM 16651	Silicone tubing, 6 × 2 mm (2 m)

## Replacement Parts

6980110	Silicone O-ring, 40 × 5 mm
6980235	Silicone O-ring, 14 × 2 mm (3/pack)
6980229	Silicone O-ring, 80 × 3 mm (2/pack)
6980225	Polycarbonate plugs (10/pack)
6980227	Polycarbonate closing cap for lid
6980236	Silicone caps (10/pack)
6980232	Polypropylene filter support (2/pack)
SM 16514E	Polycarbonate syringe filter holder, 13 mm (12/pack)

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
Sartorius AG  
Weender Landstrasse 94–108  
37075 Goettingen, Germany

Phone +49.551.308.0  
Fax +49.551.308.32 89  
www.sartorius.com

Copyright by Sartorius AG,  
Goettingen, Germany.  
All rights reserved. No part  
of this publication may  
be reprinted or translated in  
any form or by any means  
without the prior written  
permission of Sartorius AG.  
The status of the information,  
specifications and illustrations  
in this manual is indicated  
by the date given below.  
Sartorius AG reserves the  
right to make changes to the  
technology, features,  
specifications and design of  
the equipment without notice.

Status:  
April 2003, Sartorius AG,  
Goettingen, Germany

Printed in Germany on paper that has  
been bleached without any use of chlorine  
W4A000  
Publication No.: SL-6026-n03043