Fig. 1

Fig. 2

Fig. 3

Ursprungssprache

DE

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
Source language: German

EN Safety instructions

- Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Refer to the operating instructions for the drive.
- Abrasion of the milling unit or the rotating accessories can get into the medium you are working on.
- Do not use any damaged tubes e.g. with hairline cracks or splits.
- The tubes are appropriate only for a single-use. Multiple-use can destroy the tube and the milling unit.
- Excessive size and excessive hardness can damage the tubes (see **“Technical data”**).
- Twist the mill tube on securely until it clicks home audibly.
- Take care not to fill the vessel above the mark.
- Ensure that the tube is firmly attached to the bayonet lock connector on the drive unit prior to operating the device.
- After the milling operation, the tubes may be hot; please either allow it to cool down or wear appropriate personal protective equipment before detaching it from the drive.
- The processing time is limited to three minutes. The tubes must be allowed to cool down for a few minutes before processing the sample further in the same mill.
- The temperature of the material must not exceed 100 °C.
- Caution!** A sample of dry maize can reach this temperature just by milling for three minutes at maximum speed.
- The use of liquid nitrogen for directly cooling the sample in the sample vessel is **not** permitted.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Do not use the tubes in explosive atmospheres, with hazardous substances.
- Do not process materials that are radioactive, flammable or explosive, or materials that react chemically with each other releasing large amounts of energy.
- Caution!** For research purposes only! Not suitable for therapeutic or diagnostic use in human or veterinary medicine (in accordance with the IVD Directive).

Production information

- The milling tubes are for use in conjunction with the **Tube Mill control/Tube Mill 100 control**, for milling dry samples up to 5 Mohs hardness. The sample pieces must **not** be larger than 1 cm²; if the sample pieces are larger than this, they must be pre-milled before placing them in the tubes The harder the sample pieces, the smaller they must be. If for instance the hardness is 5 Mohs the edge length must be no longer than 0.5 cm (corresponding to a piece size of approx. 0.125 cm²). The tubes can be used for preparing samples for foodstuffs quality control. Commodities suitable for processing include grains, sugar, coffee and tea etc..
- Dry ice can be used to extend the application range to include moist and resilient samples such as fruits, leaves, jelly babies etc..
- The **MT 40P** may be used only for mixing powders.
- The **MTT 40/MTT 100** is intended for applications where no heavy metals (Cr or Ni) resulting from abrasion of the stainless steel beater may be permitted to contaminate the sample.
- The tubes are ideal for applications, with which no cross contamination is permissible.
- There is no need for laborious cleaning.
- The sample can be stored in the tubes.
- The tubes can be plugged and stacked using caps; these are available as accessories. Accessories see **www.ika.com**.

Application instruction

- Store the tubes in a dry place, in their packaging and at room temperature.
- Please note the limit on the usable life of sterile variants.
- Protect the tubes against ultraviolet light.
- To open the tubes, press the locking cam and twist the coupling cover anticlockwise.
- Do not fill the tubes above the maximum mark on the sample vessel. Milling less material may give better results (finer eventual particle size).
- If the samples are resilient or contain residual moisture, add a teaspoonful of dry ice (dry ice snow, not compacted dry ice) directly into the sample vessel.
- To close the tubes, twist the coupling cover clockwise until it clicks home audibly, see **Fig. 1**.
- The tube is attached by twisting it clockwise on to the bayonet coupling of the drive. The tube is correctly attached when it has clicked home audibly, see **Fig. 2**.
- When the coupling cover has been correctly screwed on to the tube and the tube is correctly attached to the drive, one of the grip cams on the sample vessel will point directly to the rear, see **Fig. 3**.
- Longer processing times will cause high temperatures in the milled material, for instance a 40 ml sample of maize milled for three minutes at maximum speed can reach a temperature of approx. 100 °C.
- If the sample is temperature-sensitive, run the mill for only a short time. Cooling the sample with dry ice snow can also prevent it overheating.
- Likewise the optimal process duration and rotating frequency must be determined by attempts.
- Before you open the tubes it is advisable to tap lightly on the mounting so that all the milled material falls into the sample vessel.
- To open the tubes, release cam and twist the coupling cover anticlockwise. Open the tube always upside down.
- Attention:** Don’t open the tube in the operation position.

Technical data

	MT 40.x	MT 100.x	MTT 40.x	MTT 100.x	MT 40P.x
Milling chamber and cover	Polypropylene (PP)				
Beater	stainless steel (AISI 301)	stainless steel (AISI 301)	Titanium	Titanium	stainless steel (AISI 301)
Seal	Polyethylene (PE)	Polyethylene (PE)	Polyethylene (PE)	Polyethylene (PE)	Polyurethane (PUR)/ Polyamide (PA)
Coupling	Polyamide (PA)				
Max. filling volume	40 ml	100 ml	40 ml	100 ml	40 ml
Max. circumferential speed	65 m/s	65 m/s	65 m/s	65 m/s	26 m/s
Max. speed	25000 rpm	25000 rpm	25000 rpm	25000 rpm	10000 rpm
Max. material size	1cm edge length (1cm²)	1cm edge length (1cm²)	1cm edge length (1cm²)	1cm edge length (1cm²)	< Ø1mm (powder)
Max. material hardness	5 Mohs				
Max. processing time	3 minutes				
Max. temperature of the medium	100 °C				
FDA conformity	Yes	Yes	Yes	Yes	No
ICH Q3D / USP chapter 232 / 233 conformity	No	No	Yes	Yes	No
Sterile variant	MT 40.x steril	MT 100.x steril	-	-	-

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