GE Healthcare Life Sciences

illustra™ ExoProStar™

Enzymatic PCR and Sequence reaction clean

Code: 29010216

Instructions for use

Warning

For research use only.

Not recommended or intended for diagnosis of disease in humans or animals.


Do not use internally or externally in humans or animals.

Recommended Protocol

- Remove the illustra Alkaline Phosphatase and Exonuclease I from the freezer and keep on ice whilst preparing the reaction.
- 2. Take a 5 µl Aliquot of the completed PCR reaction mix.
- 3. Add 1 μ l of illustra Alkaline Phosphatase to the reaction mix.
- 4. Add 1 µl of illustra Exonuclease I to the reaction mix.
- 5. Incubate at 37°C for 15 minutes.
- 6. Incubate at 80°C for 15 minutes to inactivate the enzymes.

The PCR product is now ready for use in sequencing or other down-stream applications as required.

This reaction is scalable for larger quantities of PCR reaction mix. Simply add Alkaline Phosphatase and Exonuclease I in proportion to the volume of PCR reaction mix required.



Legal

GE, imagination at work and GE Monogram are trademarks of General Electric Company.

illustra and ExoProStar and trademarks of GE Healthcare companies.

© 2011-2012 General Electric Company—All rights reserved.

First published Sept 2011

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them.

A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

For your local office contact information, visit www.gelifesciences.com/contact

GE Healthcare UK Limited Amersham Place Little Chalfont, Buckinghamshire, HP7 9NA, UK

http://www.gelifesciences.com

GE Healthcare offices:

GE Healthcare Bio-Sciences AB Björkgatan 30, 751 84 Uppsala, Sweden

GE Healthcare Europe GmbH Munzinger Strasse 5, D-79111 Freiburg, Germany

GE Healthcare Bio-Sciences Corp. 800 Centennial Avenue, P.O. Box 1327, Piscataway, NJ 08855-1327, IISA

GE Healthcare Japan Corporation Sanken Bldg. 3-25-1, Hyakunincho, Shinjuku-ku, Tokyo 169-0073, Japan

