

GE Healthcare
Life Sciences

Amersham™ Nitrocellulose Western blotting membranes


Selection guide



Amersham Nitrocellulose Western blotting membranes

GE Healthcare Life Sciences offers an extensive line of 100% pure cellulose filter paper membranes manufactured entirely from high-quality cotton linter. The membranes are manufactured and tested specifically for chromatographic and blotting applications and contain no additives that can cause contamination during the transfer step. High purity gives the membranes the excellent wicking capability and uniform capillary action that is important for obtaining clean and even transfers during blotting.

The Amersham Protran™ nitrocellulose membranes are commonly used as matrix in Western blotting applications. The membranes exhibit high affinity for proteins, exceptional blocking ability, and are compatible with a vast variety of detection methods (e.g., chemiluminescence, fluorescence, chromogenic techniques). All blotting membranes are available as sheets suitable for the most popular gel and transfer cassettes, for convenient out-of-the-box usage with minimal sheet-to-sheet variation.



Two great brands join forces


The trusted Hybond™ and Whatman™ brands combine to give Amersham Western blotting membranes, offering sustained high quality and performance.

- Optimized for chemiluminescent and fluorescent detection
- Excellent protein binding capacity over a wide size range
- Get 30% more with the new pack sizes
- Physical and chemical composition remains the same

Migration guide



New names

Whatman Protran BA79	→	Amersham Protran NC 0.1
Whatman Protran BA83	→	Amersham Protran NC 0.2
Whatman Protran BA85	→	Amersham Protran NC 0.45
Whatman Optitran BAS83	→	Amersham Protran Supported NC 0.2
Whatman Optitran BAS85	→	Amersham Protran Supported NC 0.45
Amersham Hybond ECL 0.2 μm	→	Amersham Protran Premium NC 0.2
Amersham Hybond ECL 0.45 μm	→	Amersham Protran Premium NC 0.45





Comparison of Amersham Protran NC membranes of pore sizes (A) 0.1 μm, (B) 0.2 μm, and (C) 0.45 μm.

A Western blotting model experiment for detection of ERK1/2 in a dilution series of NIH/3T3 cell lysate was created. Limit of detection is indicated by the gray arrow. In addition, the signal intensity of the captured M_r 17 000 marker protein was monitored on each membrane. A clear correlation between membrane pore size and signal intensity of the protein band could be seen.




Comparison of Amersham Protran Supported NC membranes of pore sizes (A) 0.2 μm and (B) 0.45 μm.




Membranes were subjected to the same Western blotting model experiment as described for Amersham Protran membranes. Limit of detection is indicated by the gray arrow. The signal intensity of the M_r 17 000 marker protein captured on the membranes indicates that the supported membrane with smaller poresize are best suitable for detection of proteins in the lower molecular weight range.



Comparison of Amersham Protran Premium NC membranes of pore sizes (A) 0.2 μm and (B) 0.45 μm.

Cy™3- conjugated Amersham ECL Plex™ Western Blotting system was used for detection of ERK1/2 in NIH/3T3 cell lysate. Autofluorescence profiles of Amersham Protran Premium NC 0.2, Amersham Protran Premium NC 0.45, and Amersham Hybond PVDF 0.2 membranes were monitored at the wavelengths for Cy3 and Cy5. The PVDF membrane exhibit much higher autofluorescence compared with the nitrocellulose membranes.



Product	Code number			Quantity		
	0.1 µm	0.2 µm	0.45 µm			
Amersham Protran NC 	10-6000-00	10-6000-01	10-6000-02	1 roll	300 mm × 4 m	30% more
	10-6000-05	10-6000-06	10-6000-07	1 roll	200 mm × 4 m	30% more
	10-6000-10	10-6000-11	10-6000-12	1 roll	150 mm × 4 m	30% more
	10-6000-31	10-6000-32	10-6000-33	5 sheets	300 × 600 mm	
	10-6000-45	10-6000-43	10-6000-41	25 sheets	200 × 200 mm	
	10-6000-46	10-6000-44	10-6000-42	10 sheets	200 × 200 mm	50% more
	10-6000-64	10-6000-63	10-6000-62	25 sheets	140 × 160 mm	
	10-6000-76	10-6000-74	10-6000-72	25 sheets	100 × 100 mm	
	10-6000-77	10-6000-75	10-6000-73	10 sheets	100 × 100 mm	
	10-6000-95	10-6000-94	10-6000-93	25 sheets	80 × 90 mm	
Amersham Protran Sandwich NC (Precut membrane and filter papers)	10-6001-05	10-6001-04	10-6001-03	10 + 20	140 × 160 mm	
	10-6001-16	10-6001-15	10-6001-14	10 + 20	80 × 90 mm	
Amersham Protran Supported NC 	10-6000-15	10-6000-16	10-6000-16	1 roll	300 mm × 4 m	30% more
	10-6000-17	10-6000-18	10-6000-18	1 roll	200 mm × 4 m	30% more
	10-6000-19	10-6000-20	10-6000-20	1 roll	150 mm × 4 m	
	10-6000-36	10-6000-37	10-6000-37	5 sheets	300 × 600 mm	
	10-6000-54	10-6000-52	10-6000-52	25 sheets	200 × 200 mm	
	10-6000-53	10-6000-51	10-6000-51	25 sheets	200 × 200 mm	
	10-6000-68	10-6000-67	10-6000-67	25 sheets	140 × 160 mm	
	10-6000-83	10-6000-85	10-6000-85	10 sheets	100 × 100 mm	
	10-6000-82	10-6000-84	10-6000-84	25 sheets	100 × 100 mm	
	10-6000-99	10-6000-98	10-6000-98	25 sheets	80 × 90 mm	
Amersham Protran Supported Sandwich NC (Precut membrane and filter papers)	10-6001-09	10-6001-08	10-6001-08	10 + 20	140 × 160 mm	
	10-6001-20	10-6001-19	10-6001-19	10 + 20	80 × 90 mm	
Amersham Protran Premium NC 	10-6000-04	10-6000-03	10-6000-03	1 roll	300 mm × 4 m	30% more
	10-6000-09	10-6000-08	10-6000-08	1 roll	200 mm × 4 m	30% more
	10-6000-14	10-6000-13	10-6000-13	1 roll	150 mm × 4 m	
	10-6000-35	10-6000-34	10-6000-34	5 sheets	300 × 600 mm	
	10-6000-50	10-6000-48	10-6000-48	25 sheets	200 × 200 mm	
	10-6000-49	10-6000-47	10-6000-47	10 sheets	200 × 200 mm	
	10-6000-66	10-6000-65	10-6000-65	25 sheets	140 × 160 mm	
	10-6000-81	10-6000-79	10-6000-79	25 sheets	100 × 100 mm	
	10-6000-80	10-6000-78	10-6000-78	10 sheets	100 × 100 mm	
	10-6000-97	10-6000-96	10-6000-96	25 sheets	80 × 90 mm	
Amersham Protran Premium Sandwich NC (Precut membrane and filter papers)	10-6001-07	10-6001-06	10-6001-06	10 + 20	140 × 160 mm	
	10-6001-18	10-6001-17	10-6001-17	10 + 20	80 × 90 mm	

For local office contact information, visit
www.gelifesciences.com/contact

www.gewesternblotting.com/western-blotting-product

GE Healthcare Bio-Sciences AB
 Björkgatan 30
 751 84 Uppsala
 Sweden

GE, imagination at work, and GE monogram are trademarks of General Electric Company.

Amersham, Cy, CyDye, ECL Plex, Hybond, Protran, and Whatman are trademarks of GE Healthcare companies.

CyDye™ products are manufactured under an exclusive license from Carnegie Mellon University and are covered by US patent numbers 5,569,587 and 5,627,027. The purchase of CyDye products includes a limited license to use the CyDye products for internal research and development but not for any commercial purposes. A license to use the CyDye products for commercial purposes is subject to a separate license agreement with GE Healthcare. Commercial use shall include:

1. Sale, lease, license or other transfer of the material or any material derived or produced from it.

© 2013 General Electric Company—All rights reserved.

First published Oct. 2013

All goods and services are sold subject to the terms and conditions of sale of the company within GE Healthcare which supplies them. A copy of these terms and conditions is available on request. Contact your local GE Healthcare representative for the most current information.

GE Healthcare UK Limited, Amersham Place
 Little Chalfont, Buckinghamshire, HP7 9NA, UK

GE Healthcare Europe, GmbH, Munzinger Strasse 5
 D-79111 Freiburg, Germany

GE Healthcare Bio-Sciences Corp., 800 Centennial Avenue, P.O. Box 1327
 Piscataway, NJ 08855-1327, USA

GE Healthcare Japan Corporation, Sanken Bldg., 3-25-1, Hyakunincho
 Shinjuku-ku, Tokyo 169-0073, Japan

