## **GE Healthcare**

# the smart choice for reliable results

### **Amersham ECL Direct**

Fast and convenient system for direct, nonradioactive labeling and detection. Suitable for high- to medium-target applications such as colony/plaque screens, dot blots, and PCR product analysis.



# Amersham ECL Direct

ECL Direct<sup>™</sup> is the most rapid nonradioactive labeling and detection system available. It is an easy-to-use system for high- to medium-target amounts. ECL Direct employs the convenience of direct enzyme labeling (no blocking or antibody stages) to label DNA or RNA probes with horseradish peroxidase (HRP).

### Other features include:

- Rapid results: Direct probe labeling in a 10-min reaction, 1 h from hybridization to detection with ECL Direct, Hybond<sup>™</sup> N+ and Hyperfilm<sup>™</sup> ECL.
- **Ease of handling:** Eliminates handling, waste, and regulatory issues associated with the use of radioactivity.
- Simple reprobing: No need to strip blots before reprobing.

This makes ECL Direct the ideal labeling and detection system for applications requiring rapid throughput and fast results.

# Using ECL Direct

#### **Probe labeling**

To label the probe, denatured single-stranded DNA or RNA probe is added to the labeling solution. The positively charged HRP complex is bound electrostatically to the negatively charged DNA. Glutaraldehyde is added to covalently cross-link the HRP to the probe (Fig 1). The probe can either be used directly in hybridization without purification or it can be stored in 50% glycerol for several months at -15 °C to -30 °C.




Fig 1. Probe labeling.

#### Hybridization

Specially formulated ECL Gold Buffer which stabilizes the activity of the enzyme is used for both the 15-min prehybridization and the probe hybridization step, which is typically carried out overnight (Fig 2). For colony or plaque screens, a 2-4 h hybridization is sufficient.





Fig 2. Hybridization.

#### Detection


For detection, the blots are hybridized for 1 min in the detection reagent which combines luminol—a substrate for HRP—and an enhancer to maximize light output (Fig 3). The blots are then exposed to either Hyperfilm ECL for several minutes to 1 h or a light-capture scanning device. ECL Direct detects down to 500 fg of target DNA.



# Fast and accurate results with Amersham ECL Direct

#### Strong signals with low background


ECL Direct provides consistent results with the use of ECL Gold Hybridization Buffer. The ECL Gold Hybridization Buffer contains a novel rate enhancer that gives high sensitivity with very low backgrounds on nylon membranes. It also enables lower hybridization temperatures which protect the HRP enzyme.



**Fig 4**. Random YAC DNA digested with *Pvu* II and hybridized with Alu I human DNA repetitive probe (exposure time of 30 min), using ECL Direct labeling.

#### Reprobing

HRP is inactivated by the chemiluminescent reaction so there is no need to strip blots before reprobing. This saves time and avoids possible membrane damage, which is often the limiting factor for the number of times a blot can be reprobed.



**Fig 5**. Reprobing. Comparison of 1st and 6th reprobing of Southern-blotted human genomic DNA bound to Hybond-N+. Minimal reduction in signal intensity is observed. (A) 1st reprobing; (B) 6th reprobing.

#### **Recommended applications**

ECL Direct is particularly useful for applications that require a blot to be reprobed several times and also for screening applications requiring rapid throughput and fast results. Examples include:

- High- to medium-target Southern and Northern blots
- Colony/plaque lifts
- Dot/slot blots
- Genome mapping
- PCR product analysis
- Library screens

### Ordering information

Products	Quantity	Code Number
ECL Direct Labeling and Detection System	To label 5 µg	RPN3000
ECL Direct Labeling and Detection System	To label 10 µg	RPN3001
ECL Direct Labeling module	To label 5 µg	RPN3005
ECL Detection Reagents	For 2000 cm <sup>2</sup> membrane	RPN3004
ECL Detection Reagents	For 4000 cm <sup>2</sup> membrane	RPN2105
ECL Gold Hybridization Buffer	Sufficient for 4000 cm <sup>2</sup> membrane	RPN3006
Related Products		
<b>Related Products</b> VacuGene™ XL Vacuum Blotting Unit	1	80-1266-24
	1	80-1266-24 80-6095-58
VacuGene™ XL Vacuum Blotting Unit		
VacuGene™ XL Vacuum Blotting Unit PR 648 Slot Blot Manifold, 48-well	1	
VacuGene™ XL Vacuum Blotting Unit PR 648 Slot Blot Manifold, 48-well Hybond-N+	1 see catalog for full product range	80-6095-58
VacuGene™ XL Vacuum Blotting Unit PR 648 Slot Blot Manifold, 48-well Hybond-N+ Hybridization Oven/Shaker 220/240 V	1 see catalog for full product range 1	80-6095-58 RPN2510E

General Electric Company reserves the right, subject to any regulatory approval if required, to make changes in specifications and features shown herein, or discontinue the product described at any time without notice or obligation. Contact your GE Representative for the most current information.

AlkPhos Direct, Amersham, ECF, Hybond, Hyperfilm, and Typhoon are trademarks of GE Healthcare Companies.

CDP-Star is a trademark of Tropix Inc.

ECL Direct is covered by US patent numbers 373017 and 513040 and foreign equivalents and is sold under license from Institut Pasteur.

© 2006 General Electric Company—All rights reserved. GE and GE Monogram are trademarks of General Electric Company. Asia Pacific T +85 65 62751830 F +85 65 62751829

Australasia T +61 2 8820 8299 F +61 2 8820 8200

**Austria** T 01 /57606 1613 F 01 /57606 1614

**Belgium** T 0800 73 890 F 02 416 8206

Canada T 1 800 463 5800 F 1 800 567 1008

**Central & East Europe** T +43 1 972 720 F +43 1 972 722 750

#### Denmark

T +45 70 25 24 50 F +45 45 16 2424

**Eire** T 1 800 709992 F +44 1494 542010

**Finland & Baltics** T +358 9 512 3940 F +358 9 512 39439 France T 01 69 35 67 00 F 01 69 41 98 77

**Germany** T 0800 9080 711 F 0800 9080 712

**Greater China** T +852 2100 6300 F +852 2100 6338

**Italy** T 02 26001 320 F 02 26001 399

**Japan** T 81 3 5331 9336 F 81 3 5331 9370

**Korea** T 82 2 6201 3700 F 82 2 6201 3803

Latin America T +55 11 3933 7300 F +55 11 3933 7304

Middle East & Africa T +30 210 96 00 687 F +30 210 96 00 693

Netherlands T 0800-82 82 82 1

F 0800-82 82 82 4

Norway T +47 815 65 777 F +47 815 65 666

**Portugal** T 21 417 7035 F 21 417 3184

**Russia, CIS & NIS** T +7 495 956 5177 F +7 495 956 5176

**Spain** T 902 11 72 65

F 935 94 49 65

Sweden T 018 612 1900 F 018 612 1910

**Switzerland** T 0848 8028 10 F 0848 8028 11

#### UK

T 0800 515 313 F 0800 616 927

**USA** T +1 800 526 3593

F +1 877 295 8102

www.gehealthcare.com/na\_blotting

GE Healthcare UK Limited Amersham Place Little Chalfont Buckinghamshire HP7 9NA UK

