Tech Note No. 02

# Thermo Scientific Nunc Nunclon Δ TripleFlask Culturing Technique

The Thermo Scientific Nunc Nunclon  $\Delta$  TripleFlask employs conventional flat monolayer culturing on three horizontal growth surfaces. To ensure equal distribution of cells and media in each growth chamber, prepare a homogeneous cell suspension, add to flask.

Both inoculating and cell harvesting methods are addressed in this Tech Note.

#### **Materials**

Nunc™ Nunclon™ ∆ TripleFlask (Cat. Nos. 132867 and 132913)


Medium\*: 0.2-0.4 mL/cm<sup>2</sup> culture area or 100-200 mL per TripleFlask total

Cells\*: Resuspended in complete medium

Seed at usual density (cells/cm²)

Trypsin\*: 0.1% to 0.25% ( $\pm$  10 mm EDTA) EDTA\*: 0.1-2.0% in HBSS or PBS without Ca<sup>++</sup> or Mg<sup>++</sup>

\* medium, reagents and cells as determined by previous standard culture conditions



- Prepare homogeneous cell suspension. Pour into the TripleFlask, tilting flask slightly to avoid foam or bubbles. Recommended working volume is 100-200 mL
- 2. Leave the flask in the upright position for a short time to allow equilibration of liquid in each compartment
- 2a. The flask may be canted momentarily around the connecting channel corner for facilitating the equilibration of small volumes
- 3. Quickly, but gently place the flask in the incubation position
- 4. The liquid is equally distributed over the three growth surfaces
- The flask is emptied in the same way as a conventional flask. To harvest cells add 10-15 mL Trypsin



### Method

## To Inoculate

- 1. Line up TripleFlask, standing on end with the same orientation for each flask bottom (i.e. all bottoms facing left).
  - Loosen caps, but do not remove completely.
  - Position line of flasks, leaving sufficient front work space clear. Avoid passing hands and utensils over flasks.
  - Keep caps on flasks whenever possible.
- 2. Prepare homogeneous cell suspension in a convenient vessel for dispensing.
  - Gently swirl, avoiding bubble or foam formation.
- 3. Pull first flask to be inoculated out of the line of flasks. Grip cap between little finger and palm.
  - Lift off, keeping cap interior pointed downward.
- 4. Tilt flask slightly (less than 45°), with side or bottom of flask facing palm of hand.

- 5. Pour cell suspension slowly and steadily into flask.
  - Pour cell suspension into body of flask, avoiding rim contact.
  - Medium should flow down top side of flask (surface with Nunclon imprint).
  - Avoid bubble formation.
- 6. Remount cap without touching the neck of the flask.
- 7. Stand flask on end to allow medium equilibration between compartments.
- 8. Quickly, but gently tilt flask into growth position (Nunclon imprint facing up) to distribute cells and medium equally onto each level.
- 9. Incubate as usual.

#### To Harvest

- 1. Stand flask on end. Grip flask with palm facing side or bottom
- 2. Pour medium into receptacle. Medium should flow against flask top.
- 3. Rinse monolayers with PBS or standard buffer.
  - Add rinse as you would cell inoculum.
  - Rock flask gently.
  - Drain as above.
- 4. Add 10-20 mL Trypsin or other reagent as usual. Rock to distribute evenly. Pour off excess.
- 5. Incubate at 37°C for 1-2 minutes or as usual.
- 6. Dislodge cells by tapping flask with palm of hand.
- 7. Rinse cells from flask, rocking to dislodge cells and pouring to collect harvest medium.
- 8. For other harvesting options, see Thermo Scientific Nunc Tech Note No. 3: Non-enzymatic Methods for Cell Harvesting.

Austria

+43 1 801 40 0

**Belgium** 

China

Denmark

France

Germany

India

Italy

+39 02 02 95059 or

Japan

**Netherlands** 

Nordic/Baltic countries

+358 9 329 100

**North America** +1 585-586-8800

Russia/CIS +7 (812) 703 42 15

Spain/Portugal

+34 93 223 09 18 **South America** 

**Switzerland** 

**UK/Ireland** 

**Other Asian countries** 

+852 2885 4613 **Countries not listed** 

TILSPNUNCTN02 0410