

Report No.: 8855

Task:

Application field: Biology

Material: Human hair

Feed size: 10-15 mm

Feed quantity: 500 mg

Material fibrous

specification(s):

Customer < 60 μm

requirement(s):

Subsequent DNA Desoxyribo Nuklein Acid

Solution:

analysis:

Selected Mixer Mill MM 200

instrument(s):

Configuration(s): Grinding jar MM 200 zirconium oxide 25 ml;

Grinding ball zirconium oxide ø 20 mm

Parameter(s): Frequency 25 Hz

Time: 5 min.

Achieved < 63 μm homogeneous

result(s):

Remark(s):

Recommendation: For fine grinding of different hair samples, the Mixer Mill

MM 200 is suitable under the above mentioned conditions.



Report No.: 8855

Pictures of the sample



Fig. 1: original sample



Fig. 2: after grinding in MM 200