



Task:

| Application field: | Medicine / Pharmaceuticals |
|-------------------------------|--|
| Material: | Pieces of bones |
| Feed size: | 0-5 mm |
| Feed quantity: | 300 mg – 500 mg |
| Material specification(s): | tough |
| Customer requirement(s): | Sample preparation without metalic contamination, final fineness < 100 μ m |
| Subsequent analysis: | DNA Desoxyribo Nuklein Acid |

Solution:

| Selected instrument(s): | MM 200 Mixer Mill |
|----------------------------|--|
| Configuration(s): | Grinding jar with cover 25 ml, Zirconium oxide Grinding ball \emptyset 20 mm Zirconium oxide |
| Parameter(s): | Frequency 20/s, a higher frequency can break or destroy the ceramic jars |
| Time: | 60 - 90 s |
| Achieved result(s): | predominantly < 100 µm |
| Remark(s): | The optium material quantity for a grinding jar should be approximately a third part of the jar volume. |
| Recommendation: | For fine grinding of different bone substances the Mixer Mill MM 200 is suitable under the above mentioned conditions. |