

# Data Sheet (EN)

Translation of the german original

<b>Designation, Model:</b>	<b>Order number:</b>
<b>Diaphragm pumps, chemical resistant</b> Models: <b>MPC 303 Z 230v</b>	<b>415722</b>
<b>MPC 603 E 230v</b>	<b>415721</b>

<b>Figure:</b>	<b>Description:</b>
	<p>The diaphragm pump consists of a pump body and a drive motor.</p> <p>The pump casing contains the drive unit and two pump heads. Each pump head contains a moulded diaphragm and two work valves. Both pump heads are arranged opposite each other.</p> <p>The pump heads are driven via an eccentric shaft with a connecting rod.</p>

<b>Application:</b>	<b>Materials of the medium-affecting pump parts:</b>
<p><b>The Diaphragm pumps MPC 303 Z and MPC 603 E are intended for:</b></p> <ul style="list-style-type: none"> <li>• Pumping and compressing neutral and aggressive gases and vapours according to the resistance of the indicated materials.</li> <li>• Generating a vacuum down to an ultimate pressure &lt; 5 mbar.</li> <li>• Use in physical and chemical laboratories in trade and industry.</li> <li>• Use for vacuum filtration, vacuum distillation and vacuum drying, and other vacuum technology applications.</li> </ul>	<ul style="list-style-type: none"> <li>• PTFE - Form moulded diaphragms</li> <li>• PTFE - Pump heads</li> <li>• PTFE - Vacuum hose</li> <li>• PEEK - Valves</li> <li>• EPDM - Seals</li> <li>• PVDF - Screw fittings</li> </ul>

<b>Dimensional drawing:</b>	<b>Intake pressure-/Pumping speed - Diagram:</b>												
<div style="display: flex; justify-content: space-around;"> <div style="text-align: center;"> <p>Fig. 3-1. Front view of device</p> </div> <div style="text-align: center;"> <p>Fig. 3-2. Left view of device</p> </div> </div> <p>ⓘ The main dimensions are identical for all models listed here.</p> <table border="1"> <thead> <tr> <th>Pos.</th> <th>Parameters</th> <th>Data</th> </tr> </thead> <tbody> <tr> <td>a</td> <td>Width</td> <td>≤ 250 mm</td> </tr> <tr> <td>b</td> <td>Depth</td> <td>≤ 270 mm</td> </tr> <tr> <td>c</td> <td>Height</td> <td>≤ 190 mm</td> </tr> </tbody> </table>	Pos.	Parameters	Data	a	Width	≤ 250 mm	b	Depth	≤ 270 mm	c	Height	≤ 190 mm	<p>Fig. 3-3. Pumping speed characteristics S(p)</p>
Pos.	Parameters	Data											
a	Width	≤ 250 mm											
b	Depth	≤ 270 mm											
c	Height	≤ 190 mm											

DA-415721-722

# Data Sheet (EN)

Translation of the german original

<b>Technical data:</b>			
<b>Parameter</b>	<b>MPC 603 E</b>	<b>MPC 303 Z</b>	<b>Unit</b>
<b>Pump stages</b>	1	2	-
<b>Pumping speed</b> DIN 28432	4.0 / 4.4	2.0 / 2.2	m <sup>3</sup> /h
	66 / 73	33 / 36	l / min
<b>Ultimate pressure</b>	< 55	< 5	mbar
<b>Ultimate pressure with gas ballast</b>	80	8	
<b>Max. Inlet-/Outlet pressure</b>	≤1100		mbar
<b>Intake port</b>	Small flange DN 16 KF <i>optionally enclosed</i> Hose nozzle DN 8 for hose inside diameter 8 mm		-
<b>Exhaust port</b>	Hose nozzle DN 8 for hose inside diameter 8 mm		-
<b>Ambient temperature</b>	+ 10 to + 40		°C
<b>Max. Operating gas temperature</b>	≤ + 40		
<b>Bearing</b>	maintenance-free		-
<b>Reference surface sound pressure level</b> DIN EN ISO 2151	< 44		dB (A)
<b>Voltage</b>	230		V
<b>Frequency</b>	50/60		Hz
<b>Rotation speed</b>	1500		min <sup>-1</sup>
<b>Motor power</b>	180		W
<b>Operating mode</b>	S 1		-
<b>Type of protection</b> DIN EN 60529	IP 54/40		
<b>Class of insulation</b> DIN EN 600034-1	F (160°C)		
<b>Weight net/(gross)</b>	11.2/(13.0)		kg
<b>Dimensions W/D/H</b>	250 / 270 / 190		mm

<b>Accessories:</b>			
	<b>Vacuum Control Box</b> Model: VCpro 601  for measuring and regulation of vacuum  <i>Order no.</i> <b>600100</b>		<b>Vacuum Regulator with dial gauge</b> Model: DBR-A  with hose nozzle DN 8, Connection: G1/4" - M12 x 1 to pressure adjustment - suction side  <i>Order no.</i> <b>700458</b>
	<b>Vacuum hose, Rubber</b> 18 / 8 x 5 mm  <i>Order no.</i> <b>828310-4</b>		<b>Exhaust silencer</b> Connection: 1/8"  <i>Order no.</i> <b>829901</b>

	<b>Gardner Denver Thomas GmbH</b> Am Vogelherd 20 98693 Ilmenau Germany	T +49 3677 604 0 F +49 3677 604 131 <a href="mailto:welch.emea@irco.com">welch.emea@irco.com</a> <a href="http://www.welchvacuum.com">www.welchvacuum.com</a> Customer Support +49 3677 604 0
---	--	---

We are constantly working on the further development of all our product types. Reprinting or reproduction of this manual, including extracts, is not allowed without the prior written permission of Co. Gardner Denver Thomas GmbH. All rights under the copyright laws are expressly reserved by Co. Gardner Denver Thomas GmbH. We reserve the right to make changes and amendments.

The information presented in this material is based on technical data and test results of nominal units. It is believed to be accurate and reliable and is offered as aid to help in the selection of products. It is the responsibility of the user to determine the suitability of the product for the intended use and the user assumes all risk and liability whatsoever in connection therewith. Gardner Denver Thomas GmbH does not warrant, guarantee or assume any obligation or liability in connection with this information.