

Lancer ULTIMA Series Undercounter Washers

FREE Extended Warranty (2 Year)

ULTIMA series undercounter washers offer the best labware cleaning solutions in the industry. Delivering high performance in a compact footprint. High density loading demands fewer cycles. Efficient use of water, detergents and rinsing agents minimizes the environmental impact. Energy saving construction lowers total cost of ownership.



- Optional View-In-Process window includes interior LED light to view load during the cycle
- The high performance, high temperature cleaning process is managed with preset or user defined programmable cycles
- Stainless steel top and bottom spray arms and two levels of injection washing ensure thorough spray coverage
- Dosing pump included within the price
- **FREE** 10L of Lancer clean & Lancer rinse
- **FREE** installation, commissioning and basic training
- Peace of mind service contracts available on request

**Contact
us for
a great
deal**

Code	Model	Drying	Chemical Storage	Ext. Dims (w x d x h), mm
WAS2010	810LX	NA	N	609 x 850 x 739
WAS2015	815LX	Forced- air chamber	N	609 x 850 x 739
WAS2020	820LX	Forced air chamber and injectors. HEPA filtered	Y	909 x 850 x 739

Note: Services (water, drainage and typically 3 phase power) need to be prepared prior to delivery, installation and commissioning.

Valid until 30th June 2017