

# NEW BD Falcon™ Cell Culture Multi-Flask



## Key Product Features

- ✓ Two Flask sizes available
  - 3 Layer Flask = 525cm<sup>2</sup>
  - 5 Layer Flask = 875cm<sup>2</sup>
- ✓ No special equipment required
  - Same footprint as a T175 flask saving valuable incubator space
- ✓ Access with a 10ml pipette
  - To aspirate/remove media limiting contamination risk
  - To add cells to the medium therefore less trauma to the cells
  - To remove all trypsinised cells maximising yield
- ✓ Internal Mixing port
  - Allows mixing of the cell suspension within the flask resulting in the same number of cells on each layer
  - Provides efficient equalisation of media with same volume of medium in each layer
- ✓ For production of 10<sup>7</sup> to 10<sup>9</sup> adherent cells (up to 1 gram) using a simple scale-up process starting from flasks
- ✓ Tissue Culture Treated
- ✓ Non-Pyrogenic
- ✓ Sterile (SAL of 10<sup>-6</sup>)
- ✓ DNase/RNase free
- ✓ LOT number printed on each individual flask
- ✓ Visualisation of cells on the bottom layer
- ✓ Vented cap
- ✓ Graduations (ml/layer)

Code	Description	Qty/Pack	Qty/Case	List Price
353143	3 Layer flask	2	12	£96.00
353144	5 Layer flask	1	8	£108.00

**Please contact your local Territory Sales Manager for a demonstration**