

Save
up to
57% on



SCIENTIFIC
LABORATORY
SUPPLIES

IKA® Heated Stirrers

For a limited time only SLS are offering two great deals on IKA heated stirrers. When you buy an IKA RH Basic 2 IKAMAG or RCT Basic Safety Control package you can save up to 57% on 2012 list prices!



RH Basic 2 IKAMAG stirrer

An economic magnetic stirrer with stainless steel heating plate.

Key features include:

- ✓ Soft-start motor
- ✓ Stirring quantity 10Lt H₂O
- ✓ Speed range 100 – 2,000 rpm
- ✓ Stainless steel heating plate Ø = 125 mm
- ✓ Dimensions (W x D x H) 160 x 220 x 105 mm

SLS Code	Description	2012 List Price	Offer Price
STI4452	RH Basic 2 IKAMAG stirrer	£236	£139

RCT Basic Safety Control package

The workhorse of the IKA heated magnetic stirrer range, the RCT Basic Safety Control package includes ETS-D5 contact thermometer for precise temperature control with a rod and clamp.

Key features include

- ✓ Temperature and speed setting via digital display
- ✓ Aluminium heating plate for optimal heat transfer
- ✓ Stirring quantity 20Lt H₂O
- ✓ Speed range 50 – 1,500 rpm
- ✓ Control accuracy with included ETS-D5 contact thermometer +/- 1K
- ✓ Aluminium alloy heating plate Ø 135 mm
- ✓ Dimensions (W x D x H) 160 x 270 x 85 mm

SLS Code	Description	2012 List Price	Offer Price
STI3500	RCT Basic Safety Control package	£937	£400

Save
over
£500*!



*on combined list price

With a world renowned reputation for quality, both products are available for ex stock delivery and come with the usual IKA 3 year warranty (third year provided on registration with IKA direct).

To take advantage of these offers, call your nearest Sales Office.
But hurry, stocks are limited and when they're gone, they're gone!