

## NEW Integra PipetBoy Pro – Save £113!!!

As Integra's exclusive UK distributor, SLS are very pleased to be able to launch the new PIPETBOY Pro.

Replacing the PIPETBOY Comfort, PIPETBOY Pro is a pipetting aid for use with serological pipettes. Offering all the features and benefits of the Comfort model, PIPETBOY Pro offers the added advantage of a more powerful motor for improved aspiration and an LED spotlight to ensure that the pipette is clearly lit and easily visible.



### Key features include:

- Light weight at 190g
- Rechargeable AAA NiMH batteries
- Ability to preset pipetting speed between 2-10ml/s
- Suitable for use with 2-100ml pipettes
- Filter protection to 0.45µm as standard, 0.2µm available as an option
- Options include a charging station, a combined stand for use with Vacuboy and a vacuum tweezers accessory

For a limited period only\*, the PipetBoy Pro range is available at the launch price of £199 only – a saving of £113.00 on the normal list price!!

Code	Description	Special offer price
<b>I156400</b>	PIPETBOY Pro white	£199
<b>I156401</b>	PIPETBOY Pro pink	£199
<b>I156402</b>	PIPETBOY Pro green	£199
<b>I156403</b>	PIPETBOY Pro orange	£199
<b>I156404</b>	PIPETBOY Pro blue	£199

\* Offer valid until 30th September 2010, please quote "P PRO" when you order to ensure the discount is applied. (Whilst stocks last)