



Sigma-Aldrich® Offers Compatible Vials and Closures and Simple, Easy to Use Kits*

Isn't convenience fantastic?

Compatible with a wide range of autosamplers — not only are they a great way to ensure your caps and vials fit together, they also offer great value for your money!

Order now to take advantage of these budget-busting kits.

To place an order or for further information, please contact Sigma-Aldrich Customer Service: Freephone 0800 71 71 81

Technical Service advice and applications: 0800 272 572 eurtechserv@sial.com

At least 35% cheaper than purchasing individual packs of vials, caps and septa

Cat. No.	Glass Colour	Septum	Price (£)
General use, 2mL glass vials, with aluminium crimp cap, large 6.0mm opening			
29124-U	Clear	PTFE/rubber	12.60
29125-U	Clear	PTFE/silicone	16.00
29126-U	Clear	PTFE/Silicone/PTFE	20.60
29127-U	Amber	PTFE/rubber	12.60
29128-U	Amber	PTFE/silicone	19.80
29129-U	Amber	PTFE/Silicone/PTFE	21.30

Save up to 50% by buying a kit

General use, 2mL screw thread vials (8-425, standard, 4.6mm opening), with polypropylene caps			
29104-U	Clear	PTFE/silicone	18.50
29107-U	Clear	PTFE/Silicone/PTFE	22.10
29106-U	Amber	PTFE/silicone	19.80
29108-U	Amber	PTFE/Silicone/PTFE	22.90

General use, 2mL screw thread vials (10-425, large, 6.0mm opening), with polypropylene caps			
29116-U	Clear	PTFE/rubber	25.00
29118-U	Clear	PTFE/silicone	22.90
29120-U	Clear	PTFE/Silicone/PTFE	25.00
29117-U	Amber	PTFE/rubber	26.40
29119-U	Amber	PTFE/silicone	26.40
29121-U	Amber	PTFE/Silicone/PTFE	25.00

Save around 26% when purchasing one of these packs, rather than separates

Snap Seal Vials, 2mL large (6.0mm) opening, with polypropylene cap			
29141-U	Clear	PTFE/rubber	13.80
29142-U	Clear	PTFE/silicone	19.20
29143-U	Clear	PTFE/Silicone/PTFE	20.90
29144-U	Amber	PTFE/rubber	14.30
29145-U	Amber	PTFE/silicone	19.60
29146-U	Amber	PTFE/Silicone/PTFE	21.40

*All packs of 100 vials and caps.