

## Araldite Rapid Two-Part Adhesive

Rapid two-part adhesive which is long lasting, strong, and solvent-free. Bonds to ceramics, metal, wood, chipboard, masonry, glass, dry concrete, leather, rubber, cardboard, fabric and most plastics. Resistant to water, chemical, oil, impact and temperature (-30 to 80°C). Eco-friendly packaging for easy recycling with environmentally friendly product formulation.

Code	Pack	Price
ADH1000	2 x 15mL	£17.37



## Loctite 401 Instant Adhesive

Low viscosity instant adhesive ideal for the rapid bonding of a wide range of materials, including metals, plastics and elastomers. Designed for the assembly of difficult-to-bond materials requiring uniform stress distribution and high tensile and/or shear strength. High temperature resistance from -40 to 120°C.

Code	Pack	Price
ADH1005	3g	£10.97



## SLS Lab Pro Safevac Aspiration System



Designed for laboratory waste recovery and separation of liquid and solid. Can be widely used in cell culture, DNA extraction, microplate waste removal, and any other liquid separation or recovery.

- Vacuum power can be set to different levels to achieve 1 - 15mL/s aspirating speed
- Sensitive level sensor to prevent liquid overflow during operation
- Waste bottle, tubing and handle are fully autoclavable
- Lid, tubing and handle are easy for disassembling, cleaning and autoclaving
- Wide range of adapters to fit tubes and microplates

Code	Vacuum Range, mbar	Vacuum Speed, mL/S	Volume, L	Dims, w x d x h, mm	Price
SLS4950	0 - 600	1 - 15	4	180 x 240 x 450	£1320.00



## Weigh up our range of balances!

Including portable, analytical and precision balances from OHAUS.

See pages 32-59



SCIENTIFIC  
LABORATORY  
SUPPLIES

