

Terumo® Disposable Syringes

Polypropylene body and plunger. Clear barrel for viewing of foreign materials. Individually blister packed and sterile.

Code	Capacity, mL	Description	Type	Pack	Price
SYR6200	1	Tuberculin	Concentric	100	£30.10
SYR6203	3	Luer Slip	Concentric	100	£16.38
SYR6204	5	Luer Slip	Concentric	100	£24.90
SYR6206	10	Luer Slip	Eccentric	100	£31.80
SYR6208	20	Luer Slip	Eccentric	50	£23.60
SYR6211	30	Luer Slip	Eccentric	25	£38.00
SYR6212	50	Luer Slip	Eccentric	25	£43.30



Agani™ Needles for Terumo™ Syringes

Supplied in packs of 100

Code	Gauge	Colour Code	Length, mm	Price
SYR6222	18	Pink	38	£10.15
SYR6224	19	Ivory	25	£10.01
SYR6226	19	Ivory	38	£10.02
SYR6228	19	Ivory	50	£10.02
SYR6232	20	Yellow	38	£10.60
SYR6236	21	Green	25	£10.60
SYR6238	21	Green	38	£10.60
SYR6240	21	Green	50	£10.44
SYR6242	23	Blue	16	£11.75
SYR6244	23	Blue	25	£11.25
SYR6246	23	Blue	31	£10.60
SYR6247	25	Orange	16	£11.25
SYR6248	25	Orange	25	£10.60
SYR6250	26	Brown	13	£10.60



SAMCO Glass Syringes



Glass with metal luer lock tip, fully interchangeable. Other styles available on request.

Code	Alt Ref	Capacity, mL	Price
SYR2102	G813/01	1	£7.90
SYR2104	G813/02	2	£9.20
SYR2106	G813/03	5	£12.10
SYR2108	G813/04	10	£13.40
SYR2110	G813/05	20	£16.10
SYR2114	G813/06	50	£32.80

