


Blue LO-Tox™ Teflon Coated Thermometers



Teflon coated total immersion thermometer with white back and blue filling. Teflon coating offers high levels of thermometer protection and contains any glass or filling material in the event of breakage. The FEP material is completely transparent and offers resistance to virtually all chemicals and UV radiation. Total immersion.

Code	Temp Range, °C	Graduation, °C	Length, mm	Price
THE1164	-20 to 110	1.0	305	£15.50
THE1178	-20 to 150	1.0	305	£16.61
THE1166	-10 to 210	1.0	405	£18.82




Fridge Freezer Thermometer

120mm chrome and clear front horizontal fridge or freezer thermometer with blue temperature scale on white moulding with zones indicating recommended food storage. Supplied with stainless steel clips for fixture to an appliance.

- Bright red reading column for easy reading
- No batteries required
- Heavy duty

Code	Temp Range, °C	Accuracy, °C	Dims, w x d x h, mm	Price
THE4054	-40 to 30	±1.0 (@10)	120 x 30 x 19	£5.54



Hanna Checktemp® Pocket Thermometer



Delivers high accuracy temperature measurements over a wide range without concern for breakage or condensation. Provided with Hanna's unique CAL Check™ function for accurate measurements every time. The sharp-tip probe easily penetrates semi-solid products making routine temperature checks simple and quick for both incoming and outgoing goods. It is the ideal instrument for measuring temperature according to HACCP requirements.

- Complies with EN 13485
- Save battery life with auto-off feature
- Protective probe sleeve included
- Water resistant
- Battery life approximately 2000 hours of continuous use



Code	Alt Ref	Temp Range, °C	Accuracy, °C	Resolution, °C	Dims, w x d x h, mm	Price
THE4500	HI-98501	-50 to 150	±0.2 (-30 to 120), otherwise ±0.3	0.1	50 x 21 x 185	£45.40

KEY

Throughout the catalogue the following symbols are used:



New product



Service and repairs available



Calibration service available



Made in Britain



Sustainable features



Autoclavable



No. years warranty



No. years warranty can be extended



Lifetime warranty