

Cole-Parmer® Analogue Dri-Block® Heaters Cole-Parmer®

Use these heaters as a dry, constant temperature source for various lab and research applications. All models have independent indicator lights to show main power, heater action, and over temperature cut out. Use the calibrated rotary dial to set the temperature. Dual insert analog block heater holds up to two insert blocks or one 96 well block.

- Choose a heater, then blocks to suit your application needs
- Precision engineered heaters with excellent temperature accuracy
- Versatile range of interchangeable heating blocks work with many size sample tubes and plates
- Take up less storage space – lightweight and flat top surface allow units to be stacked when not in use
- Temperature uniformity: 0.2°C at 40°C within identical blocks
- Supplied with block extraction tool



Code	Alt Ref	Model	No. Block Positions*	Temp Range, °C	Heat Output, W	Dims, w x d x h, mm	Price
BL02004	36620-12	FDB02AD	2	Ambient +5 to 100	300	202 x 260 x 105	£732.00
BL02008	36620-14	FDB03AD	3	Ambient +5 to 200	450	279 x 260 x 105	£899.00

* Single width blocks

Cole-Parmer® Digital Dri-Block® Heaters Cole-Parmer®

Use these heaters as a dry, constant temperature source for various lab and research applications. All models have independent indicator lights to show main power, heater action, and over temperature cut out. A heater that accepts two single width blocks will accept one double width block. Use three and four block models for high throughput applications such as chemical oxygen demand (COD).

- Precision engineered
- Excellent temperature accuracy
- Versatile range of interchangeable heating blocks work with many size sample tubes and plates
- Take up less storage space – lightweight with a flat top surface to allow units to be stacked when not in use
- Choice of models with 2, 3 or 4 interchangeable blocks
- Blocks sold separately



Code	Alt Ref	Model	Temp Range, °C	No. Block Positions*	Twin Control	Dims, w x d x h, mm	Price
BL01032	36620-00	DB100/2	Ambient +5 to 100	2	No	202 x 260 x 105	£739.00
BL01034	36620-02	DB100/2TC	Ambient +5 to 100	2	Yes	202 x 260 x 105	£1262.00
BL01040	36620-04	DB100/3	Ambient +5 to 100	3	No	279 x 260 x 105	£846.00
BL01042	36620-06	DB100/4	Ambient +5 to 100	4	No	356 x 260 x 105	£1130.00
BL01036	36620-08	DB200/2	Ambient +5 to 200	2	No	202 x 260 x 105	£830.00
BL01038	36620-10	DB200/3	Ambient +5 to 200	3	No	279 x 260 x 105	£887.00

* Single width blocks

KEY

Throughout the catalogue the following symbols are used:

NEW

New product

Service and repairs available

Calibration service available

Made in Britain

Sustainable features

A

Autoclavable

1

YEAR

No. years warranty

1+

YEAR

No. years warranty can be extended

L

Lifetime warranty