


Corning® Falcon® Assay Plates



- Available in a wide variety of configurations and surface chemistries
- Designed for robotisation and screening
- Optilux transparent bottom: allows direct microscopic observation of cells by standard light microscopy




96 Well, Standard Surface

- Clear

Supplied in packs of 50

Code	Well Volume, µL	Sterility	Lid	Bottom	Inner Pack	Price
353910	320	Non sterile	No	Round	5	£112.00
351172	370	Sterile	Yes	Flat	1	£163.00
351177	320	Sterile	Yes	Round	1	£163.00




96 Well, Tissue Culture Treated

- Sterile with lid
- Flat bottom
- Well volume: 370µL

Supplied in packs of 50

Code	Description	Inner Pack	Price
353296	White	5	£300.00
353072	Clear	1	£165.00
353219	Black/clear	8	£296.00
353376	Black	8	£172.00



96 Well, BioCoat™

- Sterile with lid
- Flat bottom

Supplied in packs of 50 (inner packs of 5)

Code	Description	Surface Treatment	Well Volume, µL	Price
356407	Clear	Collagen I	320	£512.00
356649	Black/clear	Collagen I	340	£850.00
356650	White/clear	Collagen I	340	£1037.00
356461	Clear	Poly-D-Lysine	370	£486.00
356640	Black/clear	Poly-D-Lysine	340	£835.00
356651	White/clear	Poly-D-Lysine	340	£833.00