




Corning® 96 Well Microplates

CORNING

- Both small volume and large volume (blocks) well formats to meet assay and storage requirements.
- Features uniform skirt heights for greater robotic gripping surface
- Solvent resistant polypropylene provides compatibility with many common organic solvents (e.g., DMSO, ethanol, methanol)
- RNase and DNase free

Code	Colour	Bottom	Well Volume, μ L	Sterility	Inner Pack	Pack	Price
3355	White	Round	360	Non sterile	25	100	£347.00
3356	Black	Round	360	Non sterile	25	100	£457.00
3359	Clear	Round	330	Sterile	25	100	£335.00
3365	Clear	Round	330	Non sterile	25	100	£298.00
3364	Clear	Flat	360	Non sterile	25	100	£298.00
3357	Clear	V	320	Sterile	25	100	£342.00
3363	Clear	V	320	Non sterile	25	100	£298.00
3343	Clear	V	450	Non sterile	10	50	£276.00
3344	Clear	V	450	Sterile	10	50	£316.00



Corning® Multiwell Plates

CORNING

- Non-reversible lids with condensation rings to reduce contamination, uniform footprint for ease in stacking
- Individual alphanumeric codes for well identification
- Treated for optimal cell attachment
- Sterilised by gamma irradiation, certified non pyrogenic
- Supplied with lid

Code	Bottom	No. of Wells	Inner Pack	Pack	Price
3506	Flat	6	5	100	£288.00
3516	Flat	6	1	50	£148.00
3512	Flat	12	5	100	£328.00
3513	Flat	12	1	50	£192.00
3524	Flat	24	1	100	£306.00
3526	Flat	24	1	50	£154.00
3527	Flat	24	5	100	£300.00
3548	Flat	48	1	100	£420.00
3596	Flat	96	1	50	£170.00
3997	Flat	96	10	50	£174.00
3598	Flat	96	5	100	£315.00
3599	Flat	96	1	100	£333.00
3585	Flat	96	5	50	£183.00
3595	Flat	96	1	50	£173.00
3696	Flat	96	1	50	£276.00
3697	Flat	96	20	100	£534.00
3799	U	96	1	50	£157.00
3894	V	96	1	50	£201.00