

Eppendorf ThermoStat™ C Block Heater/Cooler Continued

eppendorf

SmartBlocks™



Exchangeable in just seconds, these individually sensor-controlled SmartBlocks™ are the best solution for maintaining accurate, uniform temperatures. A wide variety for all common lab vessels support a broad range of applications - for use with Eppendorf ThermoMixer® C (heating/cooling and mixing) and ThermoStat™ C (heating/cooling).

- Fast and simple block exchange due to Eppendorf QuickRelease™ technology
- Optimised block design guarantees maximum temperature transfer to the sample
- Individually sensor-controlled Eppendorf SmartBlocks™ allow for block-specific calibration to provide maximum temperature accuracy and homogeneity
- Insulated thermoblocks ensure safe and ergonomic use
- Automatic block recognition for safe handling

Code	Description	Price
E5360000038	For 24 x 1.5mL reaction vessels, incl. transfer rack 1.5/2mL	£517.00
E5364000024	For 24 x reaction vessels with Ø up to 12mm, height 25-76mm	£642.00
E5366000021	For 8 x 15mL conical tubes	£528.00
E5362000035	For 24 x 2.0mL reaction vessels, incl. transfer rack 1.5/2mL	£528.00
E5309000007	For 8 x 5.0mL tubes	£517.00
E5365000028	For 4 x 50mL conical tubes	£517.00
E5307000000	For 384 well PCR plates, incl. lid	£642.00
E5306000006	For 96 well PCR plates, incl. lid	£642.00
E5318000001	For 24 x 1.5–2mL cryo tubes with max. external Ø 12.5mm, all base designs	£517.00
E5363000039	For microplates and deepwell plates, incl. lid	£642.00
E5361000031	For 24 x 0.5mL reaction vessels, incl. transfer rack 0.5mL	£517.00
E5310000002	For Eppendorf deepwell plates 96/1000µL, incl. lid	£642.00
E5316000004	For Eppendorf deepwell plates 96/500µL, incl. lid	£642.00
E3880000151	Transfer rack for 24 x 1.5/2.0 mL microtubes	£41.00
E3880000305	Transfer rack for 24 x 0.5mL microtubes	£40.50
E5363000233	Lid for SmartBlock™ 0.50–2.0mL, plates, PCR 96, PCR 384, and DWP	£116.00



SLS Select Aspirators



- Autoclavable polypropylene
- Wide shoulder handle for easy carrying and pouring
- Ideal for storage and transporting large batches of culture media, distilled water and other solutions

Code	Capacity, L	Dims, w x h, mm	Price
ASP1206	10	245 x 388	£70.00
ASP1208	20	270 x 540	£83.00