

## Lonza PC-1™ Serum-Free Medium

**LONZA**

- Low-protein medium intended for the culture of primary cells and anchorage-dependent cell lines
- Contains a complete HEPES buffering system with known amounts of insulin, transferrin, fatty acids, and proprietary proteins assembled under strict quality control procedures
- Intended for a variety of research and industrial applications
- Formulated using defined components for optimal cell growth, while maintaining the lowest possible protein content
- Complete medium, including frozen supplement, without L-glutamine



Code	Pack	Price
LZBP77232	2 x 500mL	£690.00

## Lonza Pro293™ Serum-Free Media

**LONZA**

- Optimised to support high-density growth and recombinant protein production in 293 neonatal kidney cells
- Chemically defined to ease regulatory compliance and downstream protein purification
- Very low levels of recombinant human insulin, free of animal-origin components
- Without L-glutamine or phenol red, Pluronic® F-68 is present at 0.1%

Supplied in bottles of 1L

Code	Description	Price
LZBEBP02-025Q	Pro293™S media for suspension cells	£314.00
LZBEBP12-764Q	Pro293™A media for adherent cells	£149.00



## Corning® Tissue Culture Media, Reagents and Buffers

A selection of classical media and reagents for general use in tissue culture applications. See the SLS website for full range.



Code	Description	Pack	Price
10-010-CV	Minimum Essential Media (MEM) with Earle's salts and L-glutamine	6 x 500mL	£62.00
10-013-CV	Dulbecco's Modified Eagle Medium (DMEM) with 4.5g/L glucose, L-glutamine and sodium pyruvate	6 x 500mL	£53.00
10-092-CV	Dulbecco's Modified Eagles Medium (DMEM)/Ham's F-12 (50/50 mix) with L-glutamine and 15mM HEPES	6 x 500mL	£66.00
10-040-CV	RPMI 1640 medium with L-glutamine	6 x 500mL	£53.00
10-041-CV	RPMI 1640 medium with L-glutamine and 25mM HEPES	6 x 500mL	£82.00
21-030-CV	Dulbecco's Phosphate Buffered Saline (DPBS) with calcium and magnesium	6 x 500mL	£50.40
21-031-CV	Dulbecco's Phosphate Buffered Saline (DPBS) without calcium and magnesium	6 x 500mL	£45.00
21-040-CV	Phosphate Buffered Saline (PBS) without calcium and magnesium	6 x 500mL	£46.40
21-040-CM	Phosphate Buffered Saline (PBS) without calcium and magnesium	6 x 1L	£83.00
25-051-CI	Trypsin EDTA, 0.05% Trypsin/0.53mM EDTA in HBSS	6 x 100mL	£64.00
30-001-CI	Penicillin (5000 IU) Streptomycin (5000µg/mL) 50x solution	6 x 100mL	£96.00
30-002-CI	Penicillin (10000 IU) Streptomycin (10000µg/mL) 100x solution	6 x 100mL	£105.00
46-000-CI	Molecular biology grade water tested to USP sterile purified water specifications	6 x 100mL	£62.00