

Hanna Potentiometric Titrators Continued



Code	Model	Description	Price
TIT1030	HI-931	Potentiometric titrator with 1 analogue board	£4954.00
TIT1016	HI-932-C1	Potentiometric titrator with 1 analogue board	£8255.00
TIT1018	HI-932-C2	Potentiometric titrator with 2 analogue boards	£8751.00

Accessories

Code	Alt Ref	Description	Pack	Price
TIT1020	HI-930100	Dosing pump	Each	£2374.00
TIT1022	HI-930101	Peristaltic dosing pump	Each	£2698.00
TIT1024	HI-930190	Blank burette holder	Each	£60.00
TIT1026	HI-930302	Propellers	3	£22.20
TIT1028	HI-930303	PVDF propellers	3	£51.00
TIT1040	HI-900205	5mL burette syringe	Each	£324.00
TIT1042	HI-900210	10mL burette syringe	Each	£324.00
TIT1044	HI-900225	25mL burette syringe	Each	£324.00
TIT1046	HI-900250	50mL burette syringe	Each	£283.00

2
YEARS

Hanna HI-903 Karl Fischer Volumetric Titrator



The HI-903 combines an ultra-high precision titrant delivery system with optically regulated magnetic stirring, sophisticated endpoint determination, dynamic dosing and background drift correction algorithms. The result is an extremely adaptable titrator capable of titrating with superior accuracy and precision even for samples with low moisture content. The HI-903 dispenses the titrant, detects the endpoint and performs all necessary calculations automatically.

- Powerful software and intuitive menus
- Large, colour LCD display making it simple to view results
- Equipped with a solvent handling system to reduce cell conditioning time
- Can be connected directly to your laboratory balance via serial interface
- Choose from included methods or develop a custom method for almost any application or sample type
- Using a USB flash drive or connecting the titrator to the HI 900PC application, methods (standard and user) can be upgraded, stored or deleted
- Clip-Lock™ exchangeable burette system enables users to exchange reagent burettes in a matter of seconds
- Supplied with dual platinum pin electrode, dosing pump, 5mL burette assembly with tubing, air pump assembly with tubing, beaker and bottle top assemblies and all fittings, desiccant cartridges x4 with indicating desiccant, stir bar, waste bottle, calibration key, USB cable, power cable, HI 900PC application, USB flash drive, quality certificate, ISO 8655 burette compliance report and instruction manual binder

Code	Alt Ref	Measuring Range	Resolution, ppm	Dims, w x d x h, mm	Price
TIT1032	HI-903-02	<100ppm to 100%	1 (0.0001%)	390 x 380 x 350	£5608.00