



A

Azlon® Instrument Tray



- Rectangular fabricated polypropylene tray
- Ideal for sterilising equipment
- Rigid and robust
- 22mm deep vertical sides
- Can be autoclaved at 121°C

Code	Alt Ref	Dims, l x w x d, mm	Price
TRA1054	TWR022	610 x 350 x 22	£160.00



Kartell Photographic Trays



- Deep, PVC trays with ribbed bases
- Ideal for photographic developing
- White

Code	Alt Ref	Dims, l x w x d, mm	Price
TRA3000	280	150 x 200 x 45	£14.63
TRA3002	281	200 x 250 x 60	£23.00
TRA3006	283	300 x 350 x 80	£40.80
TRA3008	284	330 x 430 x 90	£57.00
TRA3010	285	420 x 520 x 90	£81.00
TRA3014	288	420 x 540 x 180	£169.00
TRA3012	286	530 x 675 x 100	£140.00



Metal Tray Trolleys

- Constructed of tubular metal with trays and four castor tyre wheels
- Trays supplied equals the number of tiers stated
- Maximum recommended loading: 120kg
- 560mm (2 tier) or 260mm (3 tier) clearance between shelves
- Dims, h x l x w: 915 x 800 x 550mm

Code	Description	Tiers	Price
TR01000	Aluminium trays	2	£351.00
TR01002	Aluminium trays	3	£535.00
TR01004	Stainless steel trays	2	£418.00



Metal Tray Trolleys, Blue

- All steel general purpose trolleys with tubular end frames
- Fully welded
- 100mm swivel nylon castors
- Carrying capacity of 200kg in a uniformly distributed load
- Shelves have 30mm retaining lip to all 4 sides
- Supplied in sky blue (RAL 5015)
- Dimensions, l x w x h: 1060 x 460 x 950mm

Code	Tiers	Price
TR04000	2	£299.00
TR04002	3	£382.00