


## SLS Select Centrifuge Tubes



- Certified RNase and DNase free
- Sterile
- Manufactured from virgin ultra clear polypropylene
- Leak proof flat HDPE cap
- Printed graduations and a large writing area for sample identifications
- Racked tubes are packaged in recyclable card racks
- Lot traceability
- RCF: 15000xg (15mL) or 12000xg (50mL)

Supplied in packs of 500




| Code           | Capacity, mL | Bottom  | Packaging | Price   |
|----------------|--------------|---------|-----------|---------|
| <b>SLS8102</b> | 15           | Conical | Bagged    | £178.00 |
| <b>SLS8103</b> | 15           | Conical | Racked    | £199.00 |
| <b>SLS8106</b> | 50           | Conical | Bagged    | £220.00 |
| <b>SLS8110</b> | 50           | Conical | Racked    | £206.00 |
| <b>SLS8114</b> | 50           | Skirted | Bagged    | £217.00 |

## Eppendorf Centrifuge Tubes

eppendorf

- Ergonomic operation with easy open and close, and superior comfort for single-handed handling
- Stable, upright position to minimise contamination
- Superior materials and precision molding for maximal tightness and minimal sample loss
- Reliable, safe labelling due to large writing area and flat, light coloured cap
- Sterile
- Pyrogen, DNase, RNase, and DNA free (human & bacterial)
- Each individual lot is checked and certified for quality and purity
- For use up to a maximum RCF of 19500xg



| Code               | Description        | Capacity, mL | Ø, mm | Inner Pack | Pack | Price   |
|--------------------|--------------------|--------------|-------|------------|------|---------|
| <b>E0030122151</b> | Standard           | 15           | 17    | 50         | 500  | £143.00 |
| <b>E0030122178</b> | Standard           | 50           | 30    | 25         | 500  | £172.00 |
| <b>E0030122259</b> | Forensic DNA grade | 15           | 17    | 1          | 100  | £71.00  |
| <b>E0030122267</b> | Forensic DNA grade | 50           | 30    | 1          | 48   | £40.20  |

## Eppendorf 25mL Conical Tubes

eppendorf

- Reduced height gives easy sample access and requires less storage space
- Available with SnapTec® caps for single-handed opening and closing for fast liquid extraction and addition
- Screw caps have flattened and grooved sides for user-friendly ergonomic handling, stable upright positioning, and reduced contamination risk
- Centrifugation stability up to 17000xg
- Tightness is ensured from -86 up to 100°C
- Premium USP class 6 raw material and no use of slip agents, plasticizers, and biocides during manufacturing ensure highest sample integrity and optimal sample and pellet visibility
- SnapTec® cap tube can be autoclaved up to 121°C (20 min)



| Code               | Description        | Cap      | Inner Pack | Pack | Price  |
|--------------------|--------------------|----------|------------|------|--------|
| <b>E0030118405</b> | Eppendorf Quality™ | SnapTec® | 40         | 200  | £53.00 |
| <b>E0030122410</b> | Eppendorf Quality™ | Screw    | 50         | 200  | £57.00 |
| <b>E0030118413</b> | PCR clean          | SnapTec® | 40         | 200  | £57.00 |
| <b>E0030122429</b> | PCR clean          | Screw    | 50         | 200  | £61.00 |
| <b>E0030118421</b> | Sterile            | SnapTec® | 25         | 150  | £46.60 |
| <b>E0030122437</b> | Sterile            | Screw    | 25         | 200  | £64.00 |