

Corning® Large Volume Centrifuge Tubes



- Ideal for applications requiring large volume centrifugation
- Support cushions (sold separately) must be used unless the rotor has appropriately shaped V-bottom holders
- Sterile and certified non-pyrogenic
- Plug seal caps
- Maximum RCF: 6000xg

Code	Capacity, mL	Inner Pack	Pack	Price
430776	250	6	102	£251.00
431123	500	6	36	£258.00



Support Cushions

Supplied in packs of 6

Code	Description	Price
430236	250mL support cushion	£119.00
431124	500mL support cushion	£183.00

Corning® Falcon® Conical Tubes



- Easy-to-use, high quality tools to protect your precious samples
- High strength with wide temperature stability
- Meet bioanalytical-grade requirements and provide unsurpassed performance in critical research applications
- Exceptionally strong
- High clarity polypropylene or polystyrene
- Unique double threaded cap design reducing cross threading and leakage
- Easy-to-read graduations
- Consistent biological and physical properties



Code	Capacity, mL	Max RCF, xg	Material	Packaging	Pack	Price
352095	15	1800	Polystyrene	Bagged	500 (10 bags of 50)	£161.00
352099	15	1800	Polystyrene	Racked	500 (10 racks of 50)	£161.00
352096	15	12000	Polypropylene	Bagged	500 (4 bags of 125)	£140.00
352097	15	12000	Polypropylene	Racked	500 (10 racks of 50)	£175.00
352070	50	16000	Polypropylene	Bagged	500 (20 bags of 25)	£161.00
352098	50	16000	Polypropylene	Racked	500 (20 racks of 25)	£186.00
352076	175	7500	Polypropylene	Bagged	48 (4 bags of 8)	£100.00
352075	225	7500	Polypropylene	Bagged	48 (4 bags of 8)	£106.00

Support Cushions

Code	For use with	Pack	Price
352090	175 and 225mL tubes	8	£68.00