

## Polypropylene Centrifuge Tubes

- White inscription circles on caps
- Gas tight closure that cannot be overtightened
- Graduated




Code	Capacity, mL	Bottom	Dims, Ø x h, mm	Pack	Price
<b>CEN7804</b>	13	Flat	16.5 x 100	800	£358.00
<b>CEN7802</b>	14	Round	16.5 x 110	800	£380.00
<b>CEN7800</b>	15	Conical	16.5 x 120	800	£309.00
<b>CEN7808</b>	50	Self Standing	30.0 x 115	320	£161.00
<b>CEN7806</b>	50	Conical	30.0 x 115	360	£174.00

## KIMBLE® Borosilicate Glass Centrifuge Tubes



- Screw threaded tubes made from borosilicate glass
- In compliance with ASTM specification E438
- With conical shaped bottoms on the 10 and 15mL sizes and a bullet shaped bottom on the 50mL size
- Maximum RCF of 2980xg
- Supplied in a plastic shrink wrapped tray for cleanliness and safety
- Supplied without caps



Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Cap	Pack	Price
<b>TUB2700</b>	K73785/10	10	16 x 114	Screw cap, GPI thread 15-415	125	£174.00
<b>TUB2702</b>	K73785/15	15	17 x 126	Screw cap, GPI thread 15-415	125	£179.00
<b>TUB2706</b>	K73790/15	15	17 x 126	Screw cap, GPI thread 24-400	125	£176.00
<b>TUB2710</b>	K45165/15	15	17 x 118	Snap cap 2	12	£251.00
<b>TUB2704</b>	K73785/50	50	29 x 137	No cap	72	£225.00

### Caps

Code	Alt Ref	Description	For use with	Pack	Price
<b>TUB2722</b>	K73800/15415	Screw caps with PTFE-faced rubber liner, GPI thread 13-415	TUB2700, TUB2702	1000	£255.00
<b>TUB2726</b>	K73802/24400	Screw caps with PTFE-faced rubber liner, GPI thread 24-400	TUB2704	144	£226.00
<b>TUB2728</b>	K73837/2	Snap caps, polyethylene, size 2	TUB2706	500	£66.00

## SLS Select Microcentrifuge Tubes



- High transparency PP tubes for easy sample checking
- Flat cap with large frosted writing surface
- Rounded opening tab for smoother open/close routine

Supplied in packs of 1000

Code	Description	Capacity, mL	Max RCF, xg	Price
<b>CEN9304</b>	Graduated	0.5	11000	£29.00
<b>CEN9308</b>	Autoclavable	1.5	21000	£34.00

