



## CLEARLock™ Microtubes

*ClearLine*

- Polypropylene microtubes combine the security of a screw cap tube with the ease-of-use of a snap cap tube.
- Designed to withstand up to 30000xg during centrifugation without any rupture and to preserve the sample integrity
- Hinged lid prevents any unintentional opening during use
- Optimal transparency
- Temperature range: -80 to 100°C
- Autoclavable (open) at 121°C for 20 minutes
- Graduated with writing area
- Manufactured in Clean Room Class ISO 8 (UNI EN ISO 14644-1) Class 100000 (US FED STD209E)
- Human DNA, DNase, RNase, pyrogen, ATP and PCR inhibitors free certified
- CE-IVD marked

Code	Alt Ref	Volume, mL	Sterility	Pack	Price
<b>TUB1042</b>	CL05.050.0001S	0.5	Sterile	50	£32.00
<b>TUB1034</b>	CL05.001.0500	0.5	Non sterile	500	£15.00
<b>TUB1052</b>	CL05.010.1000B	0.5	Non sterile	10000	£173.00
<b>TUB1044</b>	CL15.100.0001S	1.5	Sterile	100	£59.00
<b>TUB1036</b>	CL15.002.0500	1.5	Non sterile	1000	£33.00
<b>TUB1054</b>	CL15.010.1000B	1.5	Non sterile	10000	£236.00
<b>TUB1046</b>	CL20.100.0001S	2.0	Sterile	100	£32.00
<b>TUB1038</b>	CL20.002.0500	2.0	Non sterile	1000	£36.00
<b>TUB1056</b>	CL20.008.1000B	2.0	Non sterile	8000	£213.00
<b>TUB1048</b>	CL50.050.0001S	5.0	Sterile	50	£17.00
<b>TUB1050</b>	CL50.010.0020S	5.0	Sterile	200	£45.00
<b>TUB1040</b>	CL50.002.0100	5.0	Non sterile	200	£37.00
<b>TUB1058</b>	CL50.006.0500B	5.0	Non sterile	3000	£484.00

## CLEARLine® Screw Cap Microtubes

*ClearLine*

- Manufactured from medical-grade raw material to exceed the exacting standards demanded by researchers
- Three styles of screw cap to choose from and a choice of conical bottom or self-standing tubes
- Withstand high speed centrifugation up to 16400RPM or 25000xg and temperatures from -85 to 121°C
- Autoclavable to 121°C for 20 minutes
- Tubes made in medical polypropylene and cap in medical polyethylene
- Certified DNase/RNase/pyrogen free to satisfy the highest demands in several biotechnological applications



## Sterile Tubes with O-Ring Caps

Supplied in packs of 1000 (2 bags of 500)

Code	Alt Ref	Capacity, mL	Bottom	Height, mm	Collar	Writing Area	Price
<b>TUB1060</b>	CL604/S	0.5	Conical	47.1	No	No	£119.00
<b>TUB1062</b>	CL606/S	0.5	Conical	48.5	Yes	No	£152.00
<b>TUB1064</b>	CL614/S	0.5	Self standing	47.5	No	No	£119.00
<b>TUB1066</b>	CL616/S	0.5	Self standing	49.2	Yes	No	£152.00
<b>TUB1068</b>	CL724/S	1.5	Conical	47.1	No	Yes	£129.00
<b>TUB1070</b>	CL726/S	1.5	Conical	48.2	Yes	Yes	£150.00
<b>TUB1072</b>	CL734/S	1.5	Self standing	47.6	No	Yes	£130.00
<b>TUB1074</b>	CL736/S	1.5	Self standing	49.0	Yes	Yes	£150.00
<b>TUB1076</b>	CL744/S	2.0	Conical	47.2	No	Yes	£126.00
<b>TUB1078</b>	CL746/S	2.0	Conical	48.4	Yes	Yes	£161.00
<b>TUB1080</b>	CL754/S	2.0	Self standing	48.6	No	Yes	£125.00
<b>TUB1082</b>	CL756/S	2.0	Self standing	49.8	Yes	Yes	£161.00