

SAMCO Test Tubes, Borosilicate Glass Continued



Heavy Wall, Rimless

Code	Alt Ref	Dims, Ø x h, mm	Wall Thickness, mm	Pack	Price
TES1360	G033/01	10 x 75	1.5	100	£48.20
TES1362	G033/02	12 x 75	1.5	100	£35.90
TES1364	G033/03	12 x 100	1.5	100	£41.20
TES1366	G033/05	16 x 125	1.8	100	£65.00
TES1368	G033/06	16 x 150	1.8	100	£74.00
TES1370	G033/08	18 x 150	1.8	100	£92.00
TES1372	G033/10	24 x 150	1.8	50	£79.00




SAMCO Test Tubes, Neutral Glass



- Withstand moderate temperature change
- Should not be taken from cold and placed directly into the hottest part of the flame without pre-heating
- Good chemical resistance and are suitable for use with samples which are susceptible to pH change

Medium Wall, Rimmed

Code	Alt Ref	Dims, Ø x h, mm	Wall Thickness, mm	Pack	Price
TES1204	G022/33	12 x 100	1.0	100	£27.80
TES1206	G022/34	16 x 100	1.2	100	£37.90
TES1208	G022/35	16 x 125	1.2	100	£42.70
TES1212	G022/38	18 x 150	1.2	100	£67.00
TES1214	G022/40	24 x 150	1.2	50	£38.80



Medium Wall, Rimless

Code	Alt Ref	Dims, Ø x h, mm	Wall Thickness, mm	Pack	Price
TES1220	G023/31	10 x 75	1.0	100	£22.30
TES1222	G023/32	12 x 75	1.0	100	£27.20
TES1224	G023/33	12 x 100	1.0	100	£29.40
TES1226	G023/34	16 x 100	1.2	100	£37.90
TES1228	G023/35	16 x 125	1.2	100	£42.70
TES1230	G023/36	16 x 150	1.2	100	£46.40
TES1232	G023/38	18 x 150	1.2	100	£57.00
TES1234	G023/40	24 x 150	1.2	50	£38.80

