


Kimble® Disposable Culture Tubes



- All tubes have sturdy, uniform bottoms and consistent lengths
- No marking spot
- Manufactured from soda lime glass

Supplied in packs of 1000

Code	Alt Ref	Dims, Ø x h, mm	Price
TES5100	K60AM10	10 x 75	£52.20
TES5104	K60CM13	13 x 100	£69.00
TES5108	K60EM16	16 x 100	£126.00
TES5110	K60FM165	16 x 125	£168.00
TES5112	K60GM166	16 x 150	£171.00



Thermo Scientific™ Sterilin™ Silicone Tubing



- High quality and non-toxic tubing for medical, pharmaceutical and food applications
- Aseptically manufactured
- Laser etched batch code for ease of traceability
- Ideal for use with peristaltic pumps
- Can be autoclaved or dry heat sterilised

Supplied in 15m lengths



Code	Alt Ref	Bore, mm	Wall Thickness, mm	Ext Ø, mm	Price
TUB7000	TSR0050025P	0.5	0.25	1	£83.00
TUB7002	TSR0050050P	0.5	0.5	1.5	£77.00
TUB7008	TSR0100050P	1	0.5	2	£120.00
TUB7010	TSR0100100P	1	1	3	£89.00
TUB7014	TSR0150050P	1.5	0.5	2.5	£122.00
TUB7016	TSR0150100P	1.5	1	3.5	£93.00
TUB7018	TSR0150150P	1.5	1.5	4.5	£135.00
TUB7022	TSR0200050P	2	0.5	3	£124.00
TUB7024	TSR0200100P	2	1	4	£108.00
TUB7026	TSR0200125P	2	1.25	4.5	£135.00
TUB7030	TSR0250100P	2.5	1	4.5	£105.00
TUB7032	TSR0300200P	3	2	5	£128.00
TUB7034	TSR0320080P	3.2	0.8	4.8	£137.00
TUB7036	TSR0320100P	3.2	1	5.2	£120.00
TUB7038	TSR0320160P	3.2	1.6	6.4	£157.00
TUB7042	TSR0400160P	4	1.6	7.2	£134.00
TUB7044	TSR0400200P	4	2	8	£180.00
TUB7046	TSR0450160P	4.5	1.6	7.7	£145.00
TUB7048	TSR0475160P	4.75	1.6	7.95	£168.00
TUB7050	TSR0500100P	5	1	7	£150.00
TUB7052	TSR0500160P	5	1.6	8.2	£142.00
TUB7054	TSR0550160P	5.5	1.6	8.7	£170.00
TUB7056	TSR0550250P	5.5	2.5	10.5	£235.00
TUB7058	TSR0600175P	6	1.75	9.5	£186.00
TUB7060	TSR0600200P	6	2	10	£174.00
TUB7064	TSR0630160P	6.3	1.6	9.5	£189.00
TUB7068	TSR0630320P	6.3	3.2	12.7	£328.00
TUB7070	TSR0750150P	7.5	1.5	10.5	£160.00
TUB7072	TSR0800160P	8	1.6	11.2	£207.00
TUB7074	TSR0800200P	8	2	12	£217.00