


## Thermo Scientific™ Sterilin™ Silicone Tubing Continued

**Thermo**  
SCIENTIFIC

| Code           | Alt Ref     | Bore, mm | Wall Thickness, mm | Ext Ø, mm | Price   |
|----------------|-------------|----------|--------------------|-----------|---------|
| <b>TUB7076</b> | TSR0800240P | 8        | 2.4                | 12.8      | £262.00 |
| <b>TUB7078</b> | TSR0800320P | 8        | 3.2                | 14.4      | £335.00 |
| <b>TUB7080</b> | TSR0800400P | 8        | 4                  | 16        | £497.00 |
| <b>TUB7084</b> | TSR0850400P | 8.5      | 4                  | 16.5      | £524.00 |
| <b>TUB7086</b> | TSR0900250P | 9        | 2.5                | 14        | £320.00 |
| <b>TUB7088</b> | TSR0900320P | 9        | 3.2                | 15.4      | £365.00 |
| <b>TUB7090</b> | TSR0900400P | 9        | 4                  | 17        | £453.00 |
| <b>TUB7094</b> | TSR1000200P | 10       | 2                  | 14        | £273.00 |
| <b>TUB7098</b> | TSR1270320P | 12.7     | 3.2                | 19.1      | £470.00 |
| <b>TUB7102</b> | TSR1900320P | 19       | 3.2                | 25.4      | £693.00 |



## Thermo Scientific™ Nalgene™ Non-Phthalate PVC Tubing

**Thermo**  
SCIENTIFIC

- Manufactured from the highest-grade resin materials, free from ortho-phthalate DEHP plasticisers
- Resists hardening and discoloration from a wide range of chemicals and does not impart odour or taste
- Ideal for a wide range of non-medical, general lab, research and pharmaceutical applications
- Equivalent to Tygon™ E-3603 tubing
- Marked with I.D. and regulatory compliance as well as at 10cm intervals for easy identification and measurement
- Durometer hardness is Shore A, 60
- Slips easily over glass and plastic fittings; grips securely
- Autoclavable
- Transparent

Supplied in 25m lengths

| Code           | Alt Ref   | Bore, mm | Wall Thickness, mm | Ext Ø, mm | Price   |
|----------------|-----------|----------|--------------------|-----------|---------|
| <b>TUB6268</b> | 8703-0406 | 4        | 1                  | 6         | £42.60  |
| <b>TUB6270</b> | 8703-0407 | 4        | 1.5                | 7         | £54.00  |
| <b>TUB6272</b> | 8703-0508 | 5        | 1.5                | 8         | £73.00  |
| <b>TUB6274</b> | 8703-0609 | 6        | 1.5                | 9         | £97.00  |
| <b>TUB6276</b> | 8703-0610 | 6        | 2                  | 10        | £111.00 |
| <b>TUB6278</b> | 8703-0710 | 7        | 1.5                | 10        | £116.00 |
| <b>TUB6280</b> | 8703-0811 | 8        | 1.5                | 11        | £121.00 |
| <b>TUB6282</b> | 8703-0812 | 8        | 2                  | 12        | £127.00 |
| <b>TUB6284</b> | 8703-0913 | 9        | 2                  | 13        | £134.00 |
| <b>TUB6286</b> | 8703-1014 | 10       | 2                  | 14        | £155.00 |
| <b>TUB6288</b> | 8703-1216 | 12       | 2                  | 16        | £171.00 |
| <b>TUB6290</b> | 8703-1520 | 15       | 2.5                | 20        | £76.00  |



## Gas Burner Tubing

- Neoprene rubber tubing with reinforced ends
- Ideal for use with bunsen burners

| Code           | Length, mm | Price  |
|----------------|------------|--------|
| <b>TUB2100</b> | 450        | £11.82 |
| <b>TUB2102</b> | 600        | £16.37 |
| <b>TUB2104</b> | 900        | £14.97 |