

## Tygon S3™ E-3603 Flexible Tubing Continued

| Code    | Bore, mm | Ext Ø, mm | Price   |
|---------|----------|-----------|---------|
| TUB0334 | 0.8      | 2.4       | £23.00  |
| TUB0336 | 1.6      | 3.2       | £24.00  |
| TUB0344 | 3.2      | 4.8       | £35.00  |
| TUB0346 | 3.2      | 6.4       | £134.00 |
| TUB0350 | 4.8      | 7.9       | £101.00 |
| TUB0356 | 6.4      | 11.2      | £355.00 |
| TUB0358 | 6.4      | 12.7      | £237.00 |
| TUB0364 | 7.9      | 12.7      | £203.00 |
| TUB0366 | 9.6      | 12.7      | £304.00 |
| TUB0374 | 12.7     | 17.5      | £499.00 |
| TUB0376 | 12.7     | 19.1      | £675.00 |



### Rubber Tubing

- Natural red rubber tubing
- Available with normal or heavy wall
- BS2775:1987

Supplied in 10m lengths

| Code    | Size  | Description | Wall Thickness, mm | Bore, mm | Price   |
|---------|-------|-------------|--------------------|----------|---------|
| TUB4000 | N3    | Normal wall | 0.75               | 3        | £20.07  |
| TUB4002 | N5    | Normal wall | 1.5                | 5        | £28.51  |
| TUB4004 | N6.5  | Normal wall | 1.5                | 6.5      | £32.10  |
| TUB4006 | N8    | Normal wall | 2                  | 8        | £45.56  |
| TUB4008 | N10   | Normal wall | 2                  | 10       | £54.00  |
| TUB4010 | N12.5 | Normal wall | 2.25               | 12.5     | £65.53  |
| TUB4018 | Gas   | Normal wall | 3.18               | 6.35     | £55.71  |
| TUB4050 | H3    | Heavy wall  | 3                  | 3        | £37.14  |
| TUB4052 | H5    | Heavy wall  | 4                  | 5        | £71.00  |
| TUB4054 | H6.5  | Heavy wall  | 5                  | 6.5      | £107.13 |
| TUB4056 | H8    | Heavy wall  | 6                  | 8        | £151.00 |



A

### SAMCO Transparent PVC Tubing



- Colourless, non toxic PVC tubing
- Can be steam autoclaved at 121°C
- Dimensions are to BS 2775
- Very flexible tubing which will retain its bore on bending

| Code    | Alt Ref | Bore, mm | Wall Thickness, mm | Length, m | Price  |
|---------|---------|----------|--------------------|-----------|--------|
| TUB3100 | L800/01 | 3.0      | 0.75               | 30        | £27.50 |
| TUB3102 | L800/02 | 5.0      | 1.50               | 30        | £44.70 |
| TUB3104 | L800/03 | 6.5      | 1.50               | 30        | £48.80 |
| TUB3106 | L800/04 | 8.0      | 2.00               | 30        | £72.00 |
| TUB3108 | L800/05 | 10.0     | 2.00               | 30        | £74.00 |
| TUB3110 | L800/06 | 12.5     | 2.25               | 10        | £32.20 |
| TUB3112 | L800/07 | 16.0     | 3.25               | 10        | £44.50 |
| TUB3114 | L800/08 | 20.0     | 3.25               | 10        | £42.30 |