


Kartell Unequal Tube Connectors, Straight

Kartell
LABWARE

- Inlet and outlet have different diameters
- Autoclavable polypropylene
- Tapered arms
- Dimensions are nominal, and elasticity of hoses will influence choice of size

Supplied in packs of 10

| Code | Alt Ref | Nominal OD, mm | Bore, mm | Price |
|----------------|---------|------------------|-----------|--------|
| CON7080 | 877 | 4 - 8 & 10 - 12 | 1.6 & 4.6 | £15.00 |
| CON7082 | 878 | 4 - 8 & 12 - 16 | 1.6 & 7.5 | £16.00 |
| CON7084 | 879 | 8 - 12 & 14 - 16 | 4.6 & 7.5 | £20.00 |



Kartell Disconnectors, Straight

Kartell
LABWARE

- The two pieces fit tightly yet pull apart easily
- Useful for joining hoses of glass, rubber or plastic, replacing expensive adapter
- Supplied only with the same size at each end
- With tapered and rigid arms
- Polyethylene
- Dimensions are nominal, and elasticity of hoses will influence choice of size

Supplied in packs of 10

| Code | Alt Ref | Nominal OD, mm | Bore, mm | Price |
|----------------|---------|----------------|----------|--------|
| CON7100 | 434 | 4 - 6 | 1.6 | £23.00 |
| CON7102 | 435 | 6 - 8 | 2.7 | £27.00 |
| CON7104 | 436 | 8 - 10 | 3.6 | £32.00 |
| CON7106 | 535 | 10 - 12 | 5.5 | £31.00 |
| CON7110 | 537 | 14 - 16 | 8.8 | £37.00 |



Portex® Tube Connectors, T Shape

- Rigid, straight arms
- Polypropylene

Supplied in packs of 10

| Code | Ø, mm | Price |
|----------------|-------|--------|
| CON7202 | 9 | £30.00 |
| CON7204 | 11 | £38.00 |

View our Safety & PPE catalogue

To browse our new Safety and PPE catalogue, visit the SLS website:

scientificlabs.co.uk/literature (UK & Northern Ireland)

scientificlabs.ie/literature (Republic of Ireland)

scientificlabs.com/literature (International)

