

## SAMCO Trident Vials



Trident vials are made from neutral glass under the strict conditions required by the pharmaceutical industry.



### Screw Neck Vials

Manufactured with buttress thread for optimum closure. Supplied with or without screw neck caps.

Code	Alt Ref	Form	Capacity, mL	Dims, Ø x h, mm	Neck Ø, mm	Cap	Pack	Price
<b>TUB1200</b>	T101/V1+T001/C1	Tall	1.75	11.0 x 36	5.5	Yes	882	£220.00
<b>TUB1690</b>	T101/V1	Tall	1.75	11.0 x 36	5.5	No	882	£193.00
<b>TUB1202</b>	T101/V2+T001/C1	Tall	3.50	12.5 x 46	5.5	Yes	666	£173.00
<b>TUB1692</b>	T101/V2	Tall	3.50	12.5 x 46	5.5	No	666	£153.00
<b>TUB1204</b>	T101/V3+T001/C3	Tall	7.0	17.0 x 58	9.0	Yes	399	£152.00
<b>TUB1694</b>	T101/V3	Tall	7.0	17.0 x 58	9.0	No	399	£117.00
<b>TUB1206</b>	T101/V4+T001/C3	Tall	10.5	17.0 x 67	9.0	Yes	399	£151.00
<b>TUB1696</b>	T101/V4	Tall	10.5	17.0 x 67	9.0	No	399	£113.00
<b>TUB1208</b>	T101/V5+T001/C5	Tall	14.0	20.0 x 72	13.0	Yes	264	£104.00
<b>TUB1698</b>	T101/V5	Tall	14.0	20.0 x 72	13.0	No	264	£85.00
<b>TUB1210</b>	T101/V6+T001/C5	Tall	21.25	21.0 x 86	13.0	Yes	231	£125.00
<b>TUB1700</b>	T101/V6	Tall	21.25	21.0 x 86	13.0	No	231	£81.00
<b>TUB1212</b>	T101/V7+T001/C7	Tall	28.25	23.0 x 96	15.5	Yes	190	£100.00
<b>TUB1702</b>	T101/V7	Tall	28.25	23.0 x 96	15.5	No	190	£66.00
<b>TUB1220</b>	T102/V1+T001/C5	Squat	7.0	20.0 x 42	13.0	Yes	264	£129.00
<b>TUB1704</b>	T102/V1	Squat	7.0	20.0 x 42	13.0	No	264	£71.00
<b>TUB1222</b>	T102/V2+T001/C7	Squat	14.0	23.0 x 58	15.5	Yes	190	£104.00
<b>TUB1706</b>	T102/V2	Squat	14.0	23.0 x 58	15.5	No	190	£55.00
<b>TUB1224</b>	T102/V3+T002/C3	Squat	28.25	27.5 x 72	21.0	Yes	140	£98.00
<b>TUB1708</b>	T102/V3	Squat	28.25	27.5 x 72	21.0	No	140	£51.00

### Screw Neck Caps

Black polypropylene screw caps with foil faced wads.

Code	Alt Ref	For use with	Pack	Price
<b>TUB1400</b>	T001/C1	TUB1690, TUB1692	100	£10.30
<b>TUB1402</b>	T001/C3	TUB1694, TUB1696	100	£10.30
<b>TUB1404</b>	T001/C5	TUB1698, TUB1700, TUB1704	100	£10.30
<b>TUB1406</b>	T001/C7	TUB1702, TUB1706	190	£20.10
<b>TUB1408</b>	T002/C3	TUB1708	140	£13.70