

GETINGE Lancer ULTIMA Series Labware Washing Systems Continued



ULTIMA Accessories

A large variety of standard washing accessories are available to process every type of commonly found labware. A small selection of the most popular accessories are found below. All racks used within the undercounter 810LX, 815LX & 820LX models, and freestanding models 910LX, 1300LX and 1800LXA, are interchangeable as they share a common chamber cross section. This minimises the need for additional accessories with multiple unit installations. The 1400LX/P and 1600LXP have a unique chamber footprint and as such have their own library of wash racks.

A replacement wheel kit is available to purchase for returning customers with a GETINGE Lancer washer predating the ULTIMA range. The wheel kits enable customers to re-use existing wash racks in the new ULTIMA washers. This prolongs the life cycle of accessories and helps reduce the cost when replacing a GETINGE Lancer washer.

Code	Alt Ref	Description	For use with	Price
WAS2102	PST	Basic basket, lowest level only	All models except 1400/1600	£380.00
WAS2308	GC	Reversible tray lid for PST/PSBT	All models except 1400/1600	£155.00
WAS2314	LTC	Mesh basket tall with cover (122 x 134 x 163mm)	All models except 1400/1600	£180.00
WAS2326	STC	Mesh basket short with cover (122 x 134 x 93mm)	All models except 1400/1600	£165.00
WAS2300	100 PIP	Grommet style, 100 position pipette rack	All models except 1400/1600	£1555.00
WAS2304	36 IXLC	36 jet rack, long and short, 6mm Ø	All models except 1400/1600	£1550.00
WAS2136	48 IXLC	48 jet rack, long and short, 6mm Ø	All models except 1400/1600	£1845.00
WAS2306	64 IXLC	64 jet rack, long and short, 6mm Ø	All models except 1400/1600	£2020.00
WAS2042	64 IXC	64 jet rack, short, 6mm Ø	All models except 1400/1600	£1995.00
WAS2152	GCI	Full-size injector rack cover	All models except 1400/1600	£495.00
WAS2328	36 IXLC 14	36 jet rack, short and long, 6mm Ø	1400LX/LXP	£1665.00
WAS2320	48 IXLC 14	48 jet rack, short and long, 6mm Ø	1400LX/LXP	£1980.00
WAS2318	64 IXC 14	64 jet rack, short, 6 mm Ø	1400LX/LXP	£1995.00
WAS2302	124 PIP 14	Grommet style, 124 position pipette rack	1400LX/LXP	£1680.00
WAS2036	36 IXLC 16	36 jet rack, short and long, 6mm Ø	1600LXP	£1765.00

GETINGE Lancer Cleaning Chemicals



- Developed for use in all high pressure washers
- LLL and NLL are capable of cleaning all different types of contamination including laboratory soils, "baked on" organic residues, petroleum application and other hard to clean soils
- High performance
- Contamination free - dedicated double filtration production lines
- Colour coded lids and labels for operator safety
- Phosphate free and non foaming



Supplied in packs of 10L

Code	Alt Ref	Description	Price
CLE1134	41010014	NLL neutraliser	£41.10
CLE1136	41010070	LLL detergent	£64.00

Authorised distributor



We offer a broad portfolio of products from key suppliers such as:



Honeywell

MERCK

Sigma-Aldrich®
Lab & Production Materials

Supelco®
Analytical Products

spex® an antyla scientific company