

Miele PG 8536 CD Washer Disinfectors

Miele

A glassware washer and disinfector with integrated pumps and process chemical drawer. The unit features the freely programmable Profitronic Control System and PerfectTouch control which allows easy programming by the end user, and up to 64 programme slots.

- 18 standard cleaning programmes
- DryPlus hot air drying with H13 HEPA filtration
- Spray arm sensing
- 2 integrated dispenser pumps for detergent/neutralising agent
- Drawer for 4 x 5L canisters
- Ultrasound flow control
- Optional: conductivity meter
- 3 phase 400V
- Network interface for process documentation
- Dims, w x d x h, mm: 900 x 700 x 1175



Note: services (water, drainage and 3-phase power) need to be prepared prior to delivery, installation and commissioning.

Code	Alt Ref	Description	Price
WAS1076	7831480	PG 8536 CD (for pressurised water supply)	£20653.00
WAS1078	7831490	PG 8536 CD (for unpressurised water supply)	£21188.00

Miele Baskets and Insert Accessories

Miele

Code	Alt Ref	Description	Price
WAS7200	A100	Upper basket for the connection of up to two injector modules	£520.00
WAS7202	A101	Upper basket for various inserts, height adjustable with open front vertical clearance 160mm	£331.00
WAS7204	A102	Upper basket for various inserts, height adjustable with open front vertical clearance 205mm	£383.00
WAS7206	A103	Upper basket spray arm for mesh trays and other shallow loads	£383.00
WAS7208	A150	Lower basket for the connection of up to two injector modules, automatically self-sealing docking valves	£452.00
WAS7210	A151	Open front lower basket for various DIN mesh trays and various inserts	£210.00
WAS7212	A200	Load carrier mobile unit for pipettes	£1817.00
WAS7214	A202	Two level load carrier for mesh trays	£744.00
WAS7232	A303	Module for pipettes for the optimum loading of up to 98 pipettes	£1019.00
WAS7234	A304	Module for phials for the optimum loading of up to 98 vials	£1458.00
WAS7236	A306	Module for measuring cylinders	£1790.00
A3003	A300/3	Module for optimum loading of laboratory glassware – can be equipped with 8 nozzles	£257.00
A3015	A301/5	Module for optimum loading of laboratory glassware – can be equipped with 18 nozzles	£399.00
A860	A860	Holder rack for injector modules A 301/4 and A 301/5	£142.00
A861	A861	Holder rack for injector modules A 301/4 and A 301/5	£194.00

A world of products in one place...

- Over 500,000 products
- Easy ordering capability
- Detailed product information
- Promotions and offers

