

## AstellBio Effluent Decontamination Systems

Provide your laboratory with safe, chemical free, thermal sterilisation alongside wastewater decontamination. The range consists of the AstellBio Micro EDS, which can be plumbed into existing wastewater systems; the AstellBio Sink, a wash station with automatic waste steriliser; and the AstellBio Sink and AMA440 Compact Top Loading Autoclave Combo, with both a compact autoclave and wastewater sterilising sink unit. Their automatic processing allows liquid waste to be instantly disposed of, reducing the chances of spillages and contamination, while thermal sterilisation ensures all biological material is deactivated before disposal to the drain.

- Level controlled 25L holding tank
- Level controlled 9L sterilisation tank
- Castor mounted for easily moving the units
- Optimised for low water and power consumption
- Fully programmable 5.7" colour touchscreen controller
- USB interface for downloading controller data
- Standard Astell safety features
- Full range of autoclave customisation options available for Combo option



Code	Alt Ref	Description	Price
<b>WAT1206</b>	MICR09	AstellBio Micro EDS	£35447.00
<b>WAT1204</b>	SINK9	AstellBio Sink	£35447.00
<b>WAT1202</b>	COMBO440	AstellBio Sink with Autoclave Combo	£47061.00

## Accessories

Code	Alt Ref	Description	Price
<b>WAT1208</b>	AAR130AB	Integral data printer	£1328.00
<b>WAT1210</b>	AAR122AB	Ethernet interface	£323.00
<b>WAT1212</b>	IQ/OQ001AB	IQ/OQ documents package	£1341.00
<b>WAT1214</b>	AVQ006AB	CAT III compliance	£3751.00
<b>WAT1216</b>	SOAP	Automatic soap dispenser	£1779.00

## Chemicals & Consumables



Laboratory consumables  
and chemicals to  
suit a wide range  
of applications