

## SLS Lab Pro PurA-Q+ Type 1 & 2 Water Purification Systems



Combines ultrapure and DI water in one unit saving you money and space. Type 2 water is ideal for buffer, reagent & media preparation, sample dilution and general chemistry. Type 1, dispensed direct from the unit or using an optional remote dispenser, is perfect for the most sensitive applications such as cell & tissue culture, molecular biology, electrophoresis, chromatography and spectroscopy.

- Type 2 and type 1 water quality from one unit
- Option of integrated and remote dispensers
- Intuitive, icon-based, interactive user interface with touchscreen
- Up to 2L/min ultrapure water with volumetric dispense
- Variable make up rate up to 25L/hr
- Mount on bench, under bench or on the wall
- Continuous TOC monitoring
- Low endotoxin, RNase, DNase and protease options available



Specifications	
Type 1 Water Quality, MΩ	15
Type 1 Water Rate, L/min	2
Bacteria, cfu/mL	<0.01
TOC, ppb	≤5
Feed Water Pressure, bar	0.1 to 6
Dims, w x d x h, mm	490 x 290 x 500
Weight, kg	30

Code	Model	Dispenser	Price
WAT2199	PurA-Q+ T	Integrated	£10885.00
WAT2201	PurA-Q+ R	Remote	£11352.00

### Reservoirs

Code	Capacity, L	Description	Dims, w x d x h, mm	Price
WAT2444	30	Storage reservoir	405 x 230 x 495	£932.00
WAT2446	60	Storage reservoir	480 x 330 x 495	£1089.00
WAT2168	100	Storage reservoir	480 x 360 x 810	£1588.00
WAT2190	100	Storage reservoir with pump	480 x 360 x 810	£2730.00

### Accessories and Consumables

Code	Description	Price
WAT2172	Valve kit for washing machine feed (single POU)	£163.00
WAT2452	Wall mounting kit	£108.00
WAT2448	Wall mounting kit for 30L tank	£86.00
WAT2450	Wall mounting kit for 60L tank	£105.00
WAT2174	UV upgrade kit for 30, 60 and 100L tank	£709.00
WAT2206	Remote dispenser assembly	£2131.00
WAT2178	RO membrane	£682.00
WAT2432	Sanitisation cartridge	£77.00
WAT2423	Pre-treatment cartridge - carbon	£250.00
WAT2425	Endure purification pack 1	£350.00
WAT2435	Endure purification pack 2	£534.00
WAT2430	Air vent filter for any tank	£209.00
WAT2429	POU filter (ultrafiltration)	£327.00
WAT2426	Internal UV lamp	£403.00
WAT1130	Annual consumable pack for non-ultrafiltration models - air vent filter, 2 x carbon pre-treatment cartridge, 2 x endure purification pack 1, endure purification pack 2, internal UV lamp, POU filter, sanitisation cartridge	£2240.00
WAT1132	Annual consumable pack for ultrafiltration models - air vent filter, 2 x carbon pre-treatment cartridge, 2 x endure purification pack 1, endure purification pack 2, internal UV lamp, POU filter, sanitisation cartridge	£2446.00