



## SLS Lab Pro PurA-Q Choice Type 2 or 3 Water Purification Systems



Available in two water quality versions providing general laboratory water, suitable for glass washers and autoclave feed, stability chambers, or Type 2 for buffer, reagent & media preparation, sample dilution etc. With many unique features including:

- Low cost, one box solution
- Wall or bench mounted
- Suits all general applications
- 50% recovery rate keeps water costs to a minimum
- Energy consumption <60W
- Unique disposable RO maintains peak efficiency
- Consumables changed in seconds
- Built in 35L reservoir
- Easy to install, low maintenance
- Minimal water wastage
- Mercury free UV

Specifications	
Type 2 Water Quality, MΩ-cm	1 to 15
Type 3 Water Quality, µS/cm	<40
TOC (Type 2 only), ppb	<5
Organics Rejection Rate (Type 3 only), %	<99
Feed Water Pressure, bar	0.1 to 6
Dims, w x d x h, mm	485 x 330 x 500
Weight, kg	20

Code	Model	Water Type	Make Up Rate, L/hr	Daily Usage (max), L/day	Tank UV Kit	Price
WAT2228	PurA-Q Choice 10 T2	Type 2 (DI)	10	100	No	£4340.00
WAT2230	PurA-Q Choice 10 T2 UV	Type 2 (DI)	10	100	Yes	£4333.00
WAT2220	PurA-Q Choice 10 T3	Type 3 (RO)	10	100	No	£3757.00
WAT2222	PurA-Q Choice 10 T3 UV	Type 3 (RO)	10	100	Yes	£4004.00
WAT2408	PurA-Q Choice 20 T2	Type 2 (DI)	20	200	No	£4339.00
WAT2418	PurA-Q Choice 20 T2 UV	Type 2 (DI)	20	200	Yes	£4435.00
WAT2406	PurA-Q Choice 20 T3	Type 3 (RO)	20	200	No	£4250.00
WAT2419	PurA-Q Choice 20 T3 UV	Type 3 (RO)	20	200	Yes	£4333.00

## Accessories and Consumables

Code	Description	Price
WAT2456	Wall mounting kit	£174.00
WAT2454	UV upgrade kit - tank mounted	£709.00
WAT2234	Distribution pump kit	£2214.00
WAT2172	Valve kit for washing machine feed (single POU)	£163.00
WAT2436	Pre-treatment module (RO) 20L/hr	£782.00
WAT2438	DI module	£243.00
WAT2440	Air vent filter	£225.00
WAT2442	UV lamp	£362.00
WAT1136	Annual consumable pack for WAT2228 - pre-treatment module (RO) 10L/hr, DI cartridge, air vent filter	£1091.00
WAT1138	Annual consumable pack for WAT2230 - pre-treatment module (RO) 10L/hr, 2 x DI cartridge, air vent filter, UV lamp	£1379.00
WAT1140	Annual consumable pack for WAT2220 - pre-treatment module (RO) 10L/hr, air vent filter	£704.00
WAT1142	Annual consumable pack for WAT2222 - pre-treatment module (RO) 10L/hr, air vent filter, UV lamp	£995.00
WAT1144	Annual consumable pack for WAT2408 - pre-treatment module (RO) 20L/hr, 2 x DI cartridge, air vent filter	£1133.00
WAT1146	Annual consumable pack for WAT2418 - pre-treatment module (RO) 20L/hr, 2 x DI cartridge, air vent filter, UV lamp	£1422.00
WAT1148	Annual consumable pack for WAT2406 - pre-treatment module (RO) 20L/hr, air vent filter	£746.00
WAT1150	Annual consumable pack for WAT2419 - pre-treatment module (RO) 20L/hr, air vent filter, UV lamp	£1034.00