

## Fistreem Calypso Automatic Water Stills

Offering many proven features, and available with outputs of 2 or 4L/hr. Can be bench or wall mounted. The 2L still is suitable for use in the smaller laboratory or for users with modest demands, whereas if higher outputs are required, the 4L model is ideal. Incorporating the patented 'cyclonic' vapour trap, the Calypso produces the very highest quality distilled water free from impurities including inorganics, organics, microorganisms and pyrogens. The 'cyclonic' effect created by the unique vapour trap eliminates all vapour borne droplets ensuring there is no carry over of impurities in the distillate.

- Totally enclosed powder coated cabinet with acrylic front cover provides a safe production environment
- High temperature cut-off protects the heaters in the event of water failure
- Boiler drain valve for easy cleaning without dismantling - simply introduce cleaning solution via the constant level device
- Automatically shuts off safely on a full reservoir and restarts automatically when water is drawn



Specifications	
Water, pH	5.6 - 6.0
Conductivity, $\mu\text{S}/\text{cm}$	1
Max-Min Water Pressure, kPa	70 - 560
Min Water Pressure (pre-treated), kPa	15
Dims, w x d x h, mm	260 x 325 x 790
Weight, kg	15

Code	Alt Ref	Output, L/hr	Power Rating, kW	Tap Water Flow (min), L/min	Price
<b>WAT2132</b>	WCA002	2	1.8	0.5	£2026.00
<b>WAT2136</b>	WCA004	4	3.0	1	£2034.00

### Accessories

Code	Alt Ref	Description	Price
<b>WAT2140</b>	WCA900	30L storage reservoir	£857.00
<b>WAT2142</b>	WCA901	Pre-treat feed accessory kit	£119.00

## View our Safety & PPE catalogue

To browse our new Safety and PPE catalogue, visit the SLS website:

**[scientificlabs.co.uk/literature](http://scientificlabs.co.uk/literature)** (UK & Northern Ireland)  
**[scientificlabs.ie/literature](http://scientificlabs.ie/literature)** (Republic of Ireland)  
**[scientificlabs.com/literature](http://scientificlabs.com/literature)** (International)

