



Thermo Scientific™ Sterilin™ Sampling Bottles

Thermo
SCIENTIFIC

500mL round, flexible polystyrene bottles. Available undosed or dosed with 20mg/L sodium thiosulphate, suitable for neutralising water samples with a low concentration of chlorine.

- Colour coded labels to aid differentiation samples (blue - dosed; green - undosed)
- Tamper evident cap
- Lot number for complete traceability
- Sterilised by gamma irradiation

Supplied in packs of 70

Code	Alt Ref	Description	Price
WAT4700	500WSC	Dosed	£347.00
WAT4702	500WSCNT	Undosed	£313.00



Thermo Scientific™ Orion™ AquaMate and AQUAfast

Thermo
SCIENTIFIC

A complete line of spectrophotometers, colorimeters and turbidity meters for water and wastewater analysis. Thermo Scientific™ Orion™ AquaMate spectrophotometers are the ideal instrument for water and wastewater analysis laboratories, combining innovative technology, compact design and accurate results. Portable design with uncompromised performances. Thermo Scientific™ Orion™ AQUAfast colorimeters and turbidity meters are accurate, dependable and easy to use, with handheld portability, pre-programmed methods, optimum wavelengths, user friendly displays, and easy calibration.

- Spectrophotometers that meet your requirements - choose the AquaMate 7100 Vis or AquaMate 8100 UV-Vis
- Simple and efficient sample analysis using over 260 pre-programmed methods for common reagent chemistries
- One point calibration adjustment feature ensures the highest accuracy from pre-programmed methods
- User customisable methods for absorbance, % transmittance and concentration testing
- Performance verification tests ensure wavelength accuracy and instrument functionality
- Rugged, handheld portability and waterproof options for fieldwork convenience
- User friendly interface for quick and easy results processing
- Compatibility with several tests and reagent types
- Approvals and regulatory requirements met in many industries



Code	Description	Price
AQ3010	AQUAfast IR LED turbidity meter kit	£970.00
AQ4500	AQUAfast IR & white light LED turbidity meter kit	£1604.00
AQ3140	AQUAfast COD colorimeter	£714.00
AQ3170	AQUAfast chlorine colorimeter	£523.00
AQ3700	AQUAfast colorimeter 6 wavelengths for COD & other digestion methods	£1468.00
AQ4000	AQUAfast IV colorimeter 4 wavelengths & auto ID	£1254.00
AQ7100	AquaMate 7100 Vis spectrophotometer	£4266.00
AQ8100	AquaMate 8100 UV-Vis spectrophotometer	£7405.00