



Lovibond® ET 730 Portable Flocc Tester



Idea for small mobile laboratories or field applications. It has small a footprint with 4 stirring points arranged around a central lamp for easier viewing of the flocculation process. The structure is epoxy painted and tested to give to an excellent resistance to the attack of chemical agents and corrosion.

- Non-skid base keeps the beakers stable during operation
- Integrated carrying handle
- Timer can be set for 1 - 30 minutes, or to continuous operation mode
- Illuminated centre column making it easy to visually compare samples
- Completely portable
- Height adjustable stirring rods
- Powered by a traditional wall plug, portable battery or vehicle cigarette lighter port

Code	Alt Ref	Samples	Speed Range, RPM	Dims, mm	Weight, kg	Price
WAT2064	2419150	4	20, 40, 50, 100 or 200	320 x 250 x 250	4.8	£1441.00



Lovibond® MD 600 Series Photometers



Provides the easiest and most reliable way of testing and analysing water. They provide the flexibility of a portable instrument with the quality and accuracy of a modern laboratory photometer. The MD 600 and 610 photometers are supplied with pre-programmed methods for up to 120 parameters, while the MD 640 offers additional fluorescence methods for the determination of PTSA and Fluorescein, which are used as chemical markers in water systems.

- Provide accurate, reproducible results very quickly
- Test conducted using Lovibond reagents in the form of tablets, powders, liquids, tube tests
- Fully water resistant IP68
- Large backlit graphics display
- 500 to 1000 point data storage
- Bluetooth® connectivity
- Compatible with Lovibond Aqua LX data management app

Code	Alt Ref	Model	Price
WAT2024	214020	MD 600	£1036.00
WAT2026	214025	MD 610	£1036.00
WAT2028	214140	MD 640	£1264.00



Lovibond® MD 100 Series Photometers



Photometers are portable handheld instruments for dedicated water analysis applications. A simple button layout, robust waterproof design, large backlit display and data storage function make the MD 100 easy to use, reliable and ideal for testing in the field, laboratory or plant room.

- Easy to read backlit display
- Water resistant IP68
- Compatible with liquid, tablet, powder reagents
- Supplied in a carry case with accessories and reagents
- Other parameters are available

Code	Alt Ref	Description	Price
WAT2066	276020	MD100T photometer chlorine DUO tablet	£487.00
WAT2068	278000	MD100T photometer Cl and ClO ₂	£579.00
WAT2070	276230	MD100 boiler water photometer	£547.00
WAT2072	276240	MD100 cooling water photometer	£547.00