

Xylem WTW photoLab® 7000 Spectrophotometers Continued



PC Software

Code	Alt Ref	Description	Price
WAT5111	902763	PhotoLab® Data spectral PC software	£304.00

Xylem WTW Water Test Kits



Depending on the application, there are a variety of test kits available for routine investigations. Photometers and test kits together form a system in which each is coordinated with the other depending on optics and the wavelength used.

- Quick, reliable, cost effective
- The right test for every application
- Guaranteed results through AQA/IQC
- Suitable for laboratory, infrequent work, or ease with very large sample sizes
- Round cuvette tests: reaction cuvette with barcode or method selection, 16 mm: sample addition, insertion, measurement and reading with minimum work effort, AQA support for stored results
- Reagent tests: large measurement range for 10, 20 and 50mm rectangular cuvettes, recording of the smallest concentrations in rectangular cuvettes up to 50mm, AQA support for stored results

Code	Alt Ref	Parameter	Measuring Range	Type	Pack	Price
WAT4302	250495	Ammonium (NH ₄)	0.010 to 2000mg/L NH ₄ -N, 0.01 - 2.58mg/L NH ₄ ⁺	Round cuvette test	25	£94.00
WAT4306	250424	Ammonium (NH ₄)	4.0 - 80.0mg/L NH ₄ -N, 5.2 - 103.0mg/L NH ₄ ⁺	Round cuvette test	25	£94.00
WAT4304	250329	Ammonia (NH ₃)/ Ammonium (NH ₄)	0.09 - 3.00mg/L NH ₃ (pH 8.5/25°C), 0.5 - 16.0mg/L NH ₄ -N, 0.6 - 20.6mg/L NH ₄ ⁺	Round cuvette test	25	£95.00
WAT4308	252081	Ammonia (NH ₃)/ Ammonium (NH ₄)	0.002 - 0.730mg/L NH ₃ (pH 8.5/25°C), 0.010 - 3.00mg/L NH ₄ -N, 0.013 - 3.86mg/L NH ₄ ⁺	Reagent test	250	£165.00
WAT4310	250420	Chlorine (Cl ₂)	0.03 - 6.00 Cl ₂ free & total	Round cuvette test	200*	£76.00
WAT4312	252070	Chemical Oxygen Demand (COD)	10 - 150mg/L COD (148°C, 2h)	Round cuvette test	25	£95.00
WAT4314	252071	Chemical Oxygen Demand (COD)	25 - 1500mg/L COD (148°C, 2h)	Round cuvette test	25	£95.00
WAT4316	250309	Chemical Oxygen Demand (COD)	500 - 10000mg/L COD (148°C, 2h)	Round cuvette test	25	£95.00
WAT4318	250410	Nitrate (NO ₃)	0.5 - 18.0mg/L NO ₃ -N, 2.2 - 79.7mg/L NO ₃	Round cuvette test	25	£117.00
WAT4320	252019	Nitrate (NO ₃)	23 - 225mg/L NO ₃ -N, 102 - 996mg/L NO ₃	Round cuvette test	25	£169.00
WAT4322	251993	Nitrate (NO ₃)	0.2 - 30.0mg/L NO ₃ -N, 1 - 133.0mg/L NO ₃	Round cuvette test	50	£78.00
WAT4324	252074	Nitrite (NO ₂)	0.010 - 0.700mg/L NO ₂ -N, 0.03 - 2.30mg/L NO ₂	Round cuvette test	25	£89.00
WAT4326	252096	Organic acids (volatile)	50 - 3000mg/L	Round cuvette test	25	£192.00
WAT4328	252075	Phosphate (PO ₄)	0.05 - 5.00mg/L PO ₄ -P, 0.05 - 5.0mg/L P _{ges} , 0.2 - 15.3mg/L PO ₄	Round cuvette test	25	£106.00
WAT4330	252076	Phosphate (PO ₄)	0.5 - 25.0mg/L PO ₄ -P, 0.5 - 25.0mg/L P _{ges} , 1.5 - 76.7mg/L PO ₄	Round cuvette test	25	£106.00
WAT4332	252045	Phosphate (PO ₄)	1.0 - 100.0mg/L PO ₄ -P, 3.0 - 307.0mg/L PO ₄	Reagent test	100	£151.00
WAT4334	252036	Total Organically-bound Carbon (TOC)	5.0 - 80.0mg/L TOC	Round cuvette test	25	£519.00
WAT4338	252038	Screw caps	—	—	25	£147.00

* 200 = 100 Cl₂ free + 100 Cl₂ total