

X-Ray Cassettes

- Aluminium frame, back and front plate with vinyl covering
- Gels are loaded from the front and closure affected by recessed spring-loaded catches
- Hinges and springs are made from corrosion resistant stainless steel
- Closure mechanism allows one handed operation

Code	Dims, w x l, mm	Price
MOL7036	180 x 240	£70.00
MOL7042	240 x 300	£204.00
MOL7044	300 x 400	£194.00
MOL7048	350 x 430	£129.00





X-Ray Film

- Double-sided blue-sensitive film combining fast exposure time and very good resolution
- Ideal general purpose film suitable for manual and automatic processing

Supplied in packs of 100

Code	Dims, w x l, mm	Price
MOL7016	180 x 240	£75.00
MOL7022	240 x 300	£121.00
MOL7028	350 x 430	£248.00



IT'S IN OUR DNA...



0115 982 1111



www.scientificlabs.co.uk



**SCIENTIFIC
LABORATORY
SUPPLIES**