



A

## SLS Select Wide Mouth Bottles, PP



- Autoclavable
- Polypropylene
- With cap

Code	Capacity, mL	Dims, h x Ø, mm	Pack	Price
<b>BOT2220</b>	30	62 x 35	12	£8.60
<b>BOT2222</b>	60	85 x 40	12	£11.20
<b>BOT2224</b>	125	100 x 51	12	£13.40
<b>BOT2226</b>	250	132 x 61	12	£20.80
<b>BOT2230</b>	1000	217 x 93	6	£27.60



## SLS Select Wide Mouth Bottles, LDPE



- Screw thread with cap attached
- Low density polyethylene

Supplied in packs of 6

Code	Capacity, mL	Dims, h x Ø, mm	Price
<b>SLS3548</b>	500	169 x 73	£14.90
<b>SLS3550</b>	1000	217 x 93	£23.70



## SLS Select Amber Wide Mouth Bottles, HDPE



- Amber
- High density polyethylene

Code	Capacity, mL	Dims, h x Ø, mm	Pack	Price
<b>BOT2260</b>	30	62 x 35	12	£10.20
<b>BOT2262</b>	60	85 x 40	12	£22.70
<b>BOT2264</b>	125	100 x 51	12	£16.00
<b>BOT2266</b>	250	132 x 61	12	£24.30
<b>BOT2268</b>	500	169 x 73	6	£16.10
<b>BOT2270</b>	1000	217 x 93	6	£35.50



## Thermo Scientific™ Nalgene™ Wide Mouth Lab Quality Bottles, HDPE



- Excellent chemical resistance to most acids, bases and alcohols
- Ideal for storing, shipping and packaging liquids
- High quality, laboratory grade plastic materials for dependably low leachables and extractables
- Uniform walls provide durability and resistance to splitting or punctures
- Suitable for freezer use to -100°C
- Linerless caps provide the ultimate in leakproof protection without the use of a liner that can wrinkle, cause leaks, or contaminate your reagents
- Leakproof

Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Closure Size, mm	Pack	Price
<b>BOT0080</b>	2104-0001	30	34 x 63	28	12	£12.80
<b>BOT0082</b>	2104-0002	60	39 x 86	28	12	£15.20
<b>BOT0084</b>	2104-0004	125	50 x 99	38	12	£21.00
<b>BOT0086</b>	2104-0008	250	61 x 131	43	12	£27.10
<b>BOT0088</b>	2104-0016	500	73 x 168	53	12	£44.30
<b>BOT0090</b>	2104-0032	1000	91 x 199	63	6	£29.10
<b>BOT0093</b>	2104-0048	1500	91 x 284	63	6	£40.60