

## Azlon® Wide Mouth Bottles, HDPE



- Translucent, rigid HDPE bottles with polypropylene screw caps
- Leakproof, linerless design
- Smooth inner contours for efficient pouring and ease of cleaning
- Heavy duty but lightweight and durable
- Ideal for a wide range of liquid storage and sampling applications
- Suitable for food contact applications (by EU and U.S. FDA regulations)
- Good chemical resistance

Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
<b>BOT3358</b>	BWH0030PN	30	32 x 69	10	£8.24
<b>BOT3360</b>	BWH0060PN	60	40 x 81	10	£9.90
<b>BOT3374</b>	BLH0100P	100	51 x 89	10	£10.40
<b>BOT3361</b>	BLH0125P	125	51 x 104	10	£11.90
<b>BOT3362</b>	BLH0150P	150	51 x 116	10	£12.00
<b>BOT3364</b>	BWH0250PN	250	60 x 140	10	£21.30
<b>BOT3366</b>	BWH0500PN	500	75 x 161	10	£27.20
<b>BOT3368</b>	BWH1000PN	1000	90 x 205	5	£23.70
<b>BOT3370</b>	BWH2000P	2000	120 x 242	5	£39.00



## Azlon® Graduated Wide Mouth Bottles, PP



- Clear crisp printed graduations show millilitres and fluid ounces
- Printed with large white write-on-panel for easy labelling
- Translucent polypropylene bottles with good chemical resistance
- Leakproof, linerless design
- Smooth inner contour for efficient pouring and ease of cleaning
- Suitable for food contact applications (EU and U.S. FDA regulations)
- Ideal for both liquids and solids, perfect for microbiological sampling
- Autoclavable at 121°C

Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
<b>BOT3450</b>	BWP0060PGRAD	60	40 x 82	10	£17.70
<b>BOT3452</b>	BLP0125PGRAD	125	51 x 104	10	£20.40
<b>BOT3454</b>	BWP0250PGRAD	250	60 x 141	10	£30.50
<b>BOT3456</b>	BWP0500PGRAD	500	75 x 161	10	£37.10
<b>BOT3458</b>	BWP1000PGRAD	1000	90 x 205	5	£33.80



A

## Azlon® Heavy Duty Wide Mouth Bottles, PP



- Translucent, rigid polypropylene bottles with polypropylene screw caps
- Leakproof, linerless design
- Smooth inner contours for efficient pouring and ease of cleaning
- Heavy duty and very durable in use
- Suitable for food contact applications (EU and U.S. FDA regulations)
- Ideal for a wide range of liquid storage, media preparation and sampling applications
- Good chemical resistance
- Autoclavable at 121°C

Code	Alt Ref	Capacity, mL	Dims, Ø x h, mm	Pack	Price
<b>BOT3378</b>	BWP0030PN	30	32 x 70	10	£9.30
<b>BOT3380</b>	BWP0060PN	60	40 x 82	10	£10.28
<b>BOT3372</b>	BLP0100P	100	51 x 89	10	£13.30
<b>BOT3381</b>	BLP0125P	125	51 x 104	10	£12.90
<b>BOT3376</b>	BLP0150P	150	51 x 116	10	£13.80
<b>BOT3384</b>	BWP0250P	250	140 x 60	10	£27.20
<b>BOT3386</b>	BWP0500P	500	161 x 75	10	£34.20
<b>BOT3388</b>	BWP1000P	1000	205 x 90	5	£32.50
<b>BOT3390</b>	BWP2000P	2000	242 x 120	5	£41.30



A