



### Azlon® Narrow Mouth Bottles, LDPE



- Translucent, flexible, LDPE bottles with polypropylene screw caps
- Lightweight and very durable in use
- Ideal for a wide range of liquid storage, media preparation and sampling applications
- Good general chemical resistance
- Leakproof, linerless design provides effective sealing between bottle and cap without the need for loose cap insert liners

| Code           | Alt Ref   | Capacity, mL | Dims, Ø x h, mm | Pack | Price         |
|----------------|-----------|--------------|-----------------|------|---------------|
| <b>BOT3000</b> | BNL0030PN | 30           | 32 x 70         | 10   | <b>£8.90</b>  |
| <b>BOT3002</b> | BNL0060PN | 60           | 40 x 82         | 10   | <b>£9.80</b>  |
| <b>BOT3003</b> | BNL0100PN | 100          | 51 x 90         | 10   | <b>£9.60</b>  |
| <b>BOT3004</b> | BNL0125PN | 125          | 51 x 106        | 10   | <b>£10.50</b> |
| <b>BOT3006</b> | BNL0150PN | 150          | 51 x 118        | 10   | <b>£12.40</b> |
| <b>BOT3008</b> | BNL0250P  | 250          | 60 x 141        | 10   | <b>£14.80</b> |
| <b>BOT3010</b> | BNL0500P  | 500          | 75 x 163        | 10   | <b>£20.50</b> |
| <b>BOT3012</b> | BNL1000P  | 1000         | 90 x 203        | 5    | <b>£18.90</b> |
| <b>BOT3014</b> | BNL2000P  | 2000         | 120 x 244       | 5    | <b>£34.60</b> |



### Azlon® Narrow Mouth Bottles, HDPE



- Translucent, rigid HDPE bottles with polypropylene screw caps
- Heavy duty and very durable in use
- Ideal for a wide range of liquid storage and sampling applications
- Good general chemical resistance
- Leakproof, linerless design provides effective sealing between bottle and cap without the need for loose cap insert liners

| Code           | Alt Ref   | Capacity, mL | Dims, Ø x h, mm | Closure Size, mm | Pack | Price         |
|----------------|-----------|--------------|-----------------|------------------|------|---------------|
| <b>BOT3030</b> | BNH0030PN | 30           | 32 x 70         | 18               | 10   | <b>£7.50</b>  |
| <b>BOT3032</b> | BNH0060PN | 60           | 40 x 82         | 18               | 10   | <b>£8.10</b>  |
| <b>BOT3034</b> | BNH0100PN | 100          | 51 x 90         | 23               | 10   | <b>£8.70</b>  |
| <b>BOT3036</b> | BNH0125PN | 125          | 51 x 106        | 23               | 10   | <b>£8.99</b>  |
| <b>BOT3038</b> | BNH0150PN | 150          | 51 x 118        | 23               | 10   | <b>£10.96</b> |
| <b>BOT3040</b> | BNH0250P  | 250          | 60 x 142        | 23               | 10   | <b>£13.70</b> |
| <b>BOT3042</b> | BNH0500P  | 500          | 75 x 163        | 22               | 10   | <b>£15.60</b> |
| <b>BOT3044</b> | BNH1000P  | 1000         | 90 x 203        | 22               | 5    | <b>£16.40</b> |
| <b>BOT3046</b> | BNH2000P  | 2000         | 120 x 244       | 38               | 5    | <b>£34.40</b> |



### Azlon® Narrow Mouth Bottles, PP



- High clarity, rigid polypropylene bottles with polypropylene screw caps
- Heavy duty and very durable in use
- Autoclavable at 121°C
- Ideal for a wide range of liquid storage, media preparation and sampling applications
- Good general chemical resistance
- Leakproof, linerless design provides effective sealing between bottle and cap without the need for loose cap insert liners