



## Corning® Octagonal Bottles, PET

CORNING

- PET resin makes the bottle lightweight, break resistant and as transparent as glass
- Ideal for storage and sampling of media, biologicals, and other aqueous solutions
- Packaged in two convenient, shrink wrapped trays within an outer bag to assure cleanliness
- Break resistant bottles are an ergonomic and easy to handle alternative to glass
- Bottles are sterilised to 10<sup>-6</sup> SAL and non pyrogenic
- Leakproof, tamper evident HDPE screw cap
- Molded graduations for accurate measurements
- USP class VI compliant bottles are validated non-cytotoxic and non-hemolytic

Code	Capacity, mL	Dims, w x h, mm	Graduations, mL	Pack	Price
<b>431729</b>	30	38 x 61	10	120	<b>£251.00</b>
<b>431730</b>	60	40 x 81	10	96	<b>£263.00</b>
<b>431731</b>	125	53 x 115	25	48	<b>£206.00</b>
<b>431732</b>	250	64 x 142	50	48	<b>£257.00</b>
<b>431733</b>	500	74.2 x 180	50	24	<b>£199.00</b>
<b>431734</b>	1000	96 x 220	50	24	<b>£255.00</b>



## Thermo Scientific™ Nalgene™ Heavy Duty Vacuum Bottles, PPCO

Thermo  
SCIENTIFIC

- Hold full vacuum for up to 24 hours
- Ideal for heavy-duty biotechnology lab use
- Thick walled PPCO bottles can withstand repeated room temperature vacuum applications
- Use as a waste aspirator bottle or as autoclavable vessel for scaleup activities in cell culture applications
- Large size with wide mouth
- Chemical resistant
- Rugged 53B and 83B caps with TPE gaskets

Code	Alt Ref	Capacity, L	Dims, Ø x h, mm	Closure Size	Pack	Price
<b>BOT0322</b>	2126-1000	1	93 x 226	53B	6	<b>£188.00</b>
<b>BOT0324</b>	2126-2000	2	119 x 260	53B	2	<b>£76.00</b>
<b>BOT0326</b>	2126-4000	4	155 x 338	83B	Each	<b>£62.00</b>
<b>BOT0328</b>	2126-5000	5	156 x 414	83B	Each	<b>£94.00</b>



## SLS Select Carboys

SLS

- Autoclavable polypropylene
- Wide shoulder handle for easy carrying and pouring
- Ideal for storage and transporting large batches of culture media, distilled water and other solutions

Code	Capacity, L	Colour	Dims, w x h, mm	Price
<b>ASP1200</b>	10	Translucent	249 x 390	<b>£46.10</b>
<b>ASP1204</b>	10	Amber	249 x 390	<b>£50.10</b>
<b>ASP1202</b>	20	Translucent	285 x 535	<b>£67.00</b>