



Corning® Gosselin™ Water Sampling Bottles



- Disposable blow moulded HDPE bottles
- Leakproof, tamper evident screw caps
- Available undosed without label or pre-dosed with 20mg/L sodium thiosulphate and labelled
- Break resistant and translucent
- Reverse graduations
- Sterile, Sterility Assurance Level (SAL) of 10^{-3}

Code	Alt Ref	Capacity, mL	Description	Neck Ø, mm	Pack	Price
WAT5500	FV250-05	250	Undosed	20	312	£401.20
WAT5502	HFV250-01	250	Dosed	20	312	£421.88
WAT5510	FC520-02	500	Undosed	22	100	£148.90
WAT5512	HFC520-01	500	Dosed	22	100	£156.14
WAT5515	FC1000-52	1000	Undosed	28	77	£170.61
WAT5516	HFC1000-50	1000	Dosed	28	77	£181.99



Corning® Gosselin™ Octagonal Water Sampling Bottles



- 500mL octagonal, polyethylene terephthalate (PET) sampling bottles
- Tamper evident screw caps
- Available undosed or dosed with 20mg/L sodium thiosulfate.
- Break resistant and transparent as glass
- Moulded graduations
- Label for traceability
- Sterile, Sterility Assurance Level (SAL) of 10^{-3}
- ISO 5667-3, ISO 19458
- Height with cap: 172mm
- Bottle Ø 75mm
- Internal neck Ø: 32mm

Supplied in packs of 120

Code	Alt Ref	Description	Price
WAT5551	P500B-51	Undosed	£235.76
WAT5553	HP500B-50	Dosed	£302.97



SLS Select Wash Bottles



- LDPE

Code	Description	Capacity, mL	Price
BOT9000	Jet through cap	250	£3.95
BOT9002	Jet through cap	500	£4.65
BOT9004	Jet through side	250	£2.31
BOT9006	Jet through side	400	£3.20